

WHIRL

INSTRUCTIONS FOR USE

Codes Std. Whirl 8351-8354 /
Multi Adjustable Footrest (option)



SECTION	PAGE START
1.0 STANDARD WHIRL	2
2.0 MULTI ADJUSTABLE FOOTREST (OPTION)	17



1.WHIRL

INSTRUCTIONS FOR USE

Codes 8351-8354



CONTENTS	PAGE
1.0 INTRODUCTION	3
2.0 ILLUSTRATION OF YOUR WHIRL	4
3.0 TECHNICAL DATA	4
4.0 FOR YOUR SAFETY	5
5.0 SETTING UP YOUR STANDER	6
6.0 ADJUSTING YOUR STANDER	6
7.0 TRANSFER	11
8.0 CARE AND MAINTENANCE	10
9.0 WARRANTY AND SERVICE	13
10.0 CONTINUOUS IMPROVEMENT	13
11.0 SERVICE INSPECTION	14

1.0 INTRODUCTION

Thank you for choosing your new Smirthwaite Whirl Stander.

The Whirl has been designed to give excellent symmetry and enable children from a young age to gradually develop and maintain their ability to weight bear. It is suitable for children who require minimal to high levels of support and control.



IMPORTANT!

These instructions should be read by all therapists and carers using the equipment and should be retained for future reference. The product should always be used under adult supervision.

Any incorrect use of the product and failure to follow the instructions may put the user at risk or impede the function. If you have any queries using this product or wish for further copies, please do not hesitate to contact Customer Service department on T: +44 (0) 1626 835552.


This product has been designed and manufactured specifically and solely for use by clients with special needs.



2.0 ILLUSTRATION OF YOUR WHIRL STANDER

- A. Padded thoracic support
- B. Padded pelvic support
- C. Activity tray and bowl
- D. Knee cups
- E. Mechanical angle adjustor*
- F. Indexable footrests

*An electric actuator version is also available.



3.0 TECHNICAL DATA

Size	User Height (mm)*	Max Child Weight (kg)
1	630-1020	70
2	760-1270	90

* Measured from footplate to Thoracic Support (top)

Size	1	2
Code	8351 ¹ /8353	8352 ¹ /8354
User weight limit (kg)	70	90
Thoracic support (top) to foot plate height	*630-1020	*760-1270
Pelvic support (top) to foot plate height	*410-830	*560-1030
Knee support (top) to foot plate height	*130-420	*165-520
Foot plate to floor height (excluding sandals)	70	70
Max width (between thoracic supports)	300	460
Tray size (depth and width)	610 x 690	700 x 800
Angle (°)	0 to 30	0 to 30
Base frame	890 x 660	1020 x 710
Shipping weight (including packaging) (kg)	38	40



4.0 FOR YOUR SAFETY

	STOP!
	Please read these instructions CAREFULLY and THOROUGHLY

- The user should **NOT** be left unattended whilst using the Whirl. Always ensure a responsible adult or carer is in attendance.
- Regular maintenance checks and cleaning are essential for the safe use of this equipment (see care and maintenance section).
- Use only on a flat surface and do not exceed the axilla heights or weight limits specified above.
- If any part is loose, damaged or functioning incorrectly, do not use until rectified.
- Check that all 4 castor brakes are locked and all parts are secure before the user is transferred to the stander.
- The Whirl is supplied fully assembled except for the tray. Check before use that all parts are present and undamaged. Please read the instructions below carefully to ensure safe usage.

	STOP!
	LAP STRAPS & HARNESSSES SAFETY NOTICE Lap straps and harnesses must be appropriate and safe for the user and the users clothing. Lap straps and harnesses must be checked every time the chair is used to ensure they are fitted as prescribed by the clinician, take account of the users clothing and are tightened so that the user cannot sustain injury. Checking the fit of lap straps and harnesses must be done with the user in the chair and should be undertaken as soon as the users sits in the chair.

	STOP!
	If you believe this product to be faulty – DO NOT USE – Contact Smirthwaite Ltd on T: +44 (0) 1626 835552



5.0 SETTING UP YOUR STANDER



STOP!

The user should **NOT** be left unattended whilst using the product. Always ensure a responsible therapist or carer is in attendance.

- Ensure all 4 castors are brake locked.
- Remove the activity tray
- Take the intended user's measurements.
- Set the stander up approximately to these measurements and make final adjustments when the user is in place.

Smirthwaite offer, at no extra charge, an appointment with a product specialist who will assist in the initial product set up.

6.0 ADJUSTING YOUR STANDER




IMPORTANT!

The user should **NOT** be left unattended whilst using the product. Always ensure a responsible therapist or carer is in attendance.

Provided with your Whirl you will find a 4mm and 6mm Hexagon key. This is for use to help adjust your stander.

6.1 MAIN COLUMN HEIGHT

- The main column height can be adjusted by loosening screw (A) with the 6mm hexagon key provided.
- Slide the inner column to the desired height and refasten (A) securely.
- The inner column is prevented from being removed completely by a mechanical stop mechanism.



6.2 KNEE BLOCK ADJUSTMENT

Depth (Flexion / Extension)


- Loosen the lever lock (B)
- Slide the knee block through the clamp to the desired depth.
- Retighten the lever lock (B)

Height

- Loosen screw (C) with the 6mm hexagon key
- Slide the clamp to the desired height
- Retighten the screw (C)

Lateral (Adduction/Abduction) Position

- Loosen the lever lock (D)
- Slide the knee block to the desired position, ensuring maximum contact is made
- Retighten the lever lock (D)
- Repeat for the second knee block
- By loosening the lever lock (D), it is also possible to swivel the knee block to the desired position.
- Repeat for second knee block



6.3 PELVIC SUPPORT

Height

- Loosen the 6mm screw (G)
- Slide the pelvic support section to desired height
- Retighten the 6mm screw (G)

Depth


- Loosen both lever locks (H)
- Adjust the depth of the pelvic support to the desired position
- Retighten both level locks (H)

Rotation

- Rotation of the pelvic support can be achieved by independently loosening lever locks (H) and adjusting each to suit.
- Retightened lever locks (H) when finished

Width

- Loosen 4 screws (I) using the hexagonal driver
- Slide the aluminum plates to the desired position
- Retighten 4 screws (I) to lock into position



6.4 THORACIC SUPPORT

Width


- Loosen 4 screws (K) using the hexagonal driver
- Slide the aluminum sides to the desired width
- Retighten the 4 screws (K)

Fine Height

- The wooden block (L) on the rear of the thoracic support has two screws that, when loosened, permit the movement of the thoracic support through a fine range upwards and downwards.
- Ensure the screws are fastened tight once the desired height is achieved.

Angle

- Rotate hand wheel (M) to give the desired angle.




6.5 TRAY

Height

- Loosen lever lock (N)
- Lift tray to the desired height
- Retighten level lock (N)
- Loosen 2 screws (P) using the hexagonal driver. This permits further tray height adjustment.
- Ensure screws (P) are secured when desired height is set.

Release


- To release the tray from the main column assembly, pull both plunger pins (O) and lift tray.
- When replacing the tray into position, ensure the plunger pins are securing engaged



6.6 COLUMN ANGLE ADJUSTMENT (MECHANICAL)

Angle Adjustment


- Ensure hand wheel (Q) is tight to engage the drive
- Rotate handle (R) until desired angle is achieved
- Disengage hand wheel (Q) to ensure the angle adjustment is secure



6.7 COLUMN ANGLE ADJUSTMENT (ELECTRIC)

Angle Adjustment

- The electric Whirl variant is fitted with a battery pack, actuator (S) and handset/controller (T).
- Press the Up or Down arrow on the handset to adjust stander to the desired angle.




6.8 BATTERY AND SAFETY CUT OUT

Battery Care

- Charge the battery 24 hours before use
- Keep the battery fully charged
- Charge once a week, overnight for instance
- To charge, plug the charging unit into the mains - ensure the red stop button is released
- Connect its lead to the handset (U)
- Ensure the handset is connector to the control box/battery
- To cut the power in case of emergency, simply push the red stop button (V)
- Ensure the red emergency stop (V) is twisted out

Note:


- The red stop button is used in case of emergency; do not use as an on/off switch
- Twist the red button clockwise to reset
- If the battery becomes discharged to a critical level, a warning tone will sound. Charge the battery as soon as this tone is heard
- Batteries should be replaced at least every 4 years.



6.9 INDEXING FOOT RESTS

Angle Adjustment

- Accurate foot placement can be achieved by making incremental adjustments of **each** foot rest
- Loosen the locking hand wheel (W) on the foot rest
- Pull the footrest outward slightly and rotate to the desired angle, using the stopper pin as a guide (X)
- When the desired position is achieved, simple lock the footrest securely by retightening the hand wheel (W)



7.0 TRANSFER

- Open all the straps on the Thoracic Support, Pelvic Support and Sandals
- Unlock the castors and move the Whirl in front of the seated child
- Lock the 4 castors at 45° to each other
- Release the catch located under each footrest and lift upwards, towards the child's feet
- Secure the child's feet in the sandals, using the straps provided
- The child may then be carefully positioned vertically into the stander
- Support the child, whilst securing the straps on the Thoracic Support and Pelvic Support
- Make adjustments to the Thoracic, Pelvic and Knee support as needed
- Follow the adjustment instructions as needed.



8.0 CARE AND MAINTENANCE

	IMPORTANT!
	CLEANING IS RECOMMENDED ON A REGULAR BASIS

THE CLIENT SHOULD NOT BE IN THE STANDER WHILE THE CHECKS ARE CARRIED OUT.

- Clean upholstery and wooden parts with a damp cloth and mild detergent
- Make sure brakes are released before moving
- Store the stander in a cool dry place out of direct sunlight.
- Anti-bacterial wipes may be used to clean the Whirl Tray.
- For further information please refer to MHRA or NHS cleaning guidelines

8.1 DAILY & MONTHLY CHECKS

- Check the product for damage or loose components.
- Check the structure of the product for signs of damage wear and tear.
- Keep all parts clean.

8.2 SERVICE INTERVAL

The Whirl should be serviced every 1 year. Servicing must only be undertaken by a Smirthwaite service engineer, or by a Smirthwaite trained representative.

8.3 NOMINAL SERVICE LIFE

Your product has a nominal service life of 7 years, during which full post-sales support will be available with regard to spares and servicing.

Product service life has been determined based upon the design complexity of product, and the anticipated exposure to normal use.

Good practice dictates all Smirthwaite products have been designed and manufactured to high levels of safety and quality, and will meet requirements of normal use when maintained in line with our servicing recommendations.





STOP!

If the product has been out of use for an extended period of time (6 months or more) it should always be serviced prior to being reissued.

If the product has been subjected to 'heavy' or 'constant' use, the service should be reduced to half the recommended period.

Constant and/or heavy use is considered to be:

- Daily use above 7 hours duration
- Weekly use above 5 days duration
- Monthly use above 10 months per year
- Use by a client who is at 90% to 100% of the maximum weight limit of the product. The maximum weight limit must **NEVER** be exceeded
- Use by a client who is extremely active, either voluntarily or involuntarily

8.4 EXTENDING NOMINAL SERVICE LIFE

At Smirthwaite we are proud to produce products that have a reputation for quality and durability.

We believe our products have the potential to provide benefits to our clients beyond the nominal service life documented above.

We will continue to provide full support beyond the nominal service life provided the following conditions are met:

1. A full service schedule has been maintained.
2. A full service and inspection is undertaken at the end of the nominal service life period
3. The product is subsequently serviced annually (or biannually if under 'heavy/constant' use conditions)
4. Smirthwaite reserve the right to limit support where parts/components are no longer available

8.5 DOCUMENTATION/RECORDS

- It is the responsibility of the current equipment owner to ensure the 'Instructions for Use' manual and any further manuals for accessories fitted to the equipment are handed over to the new owner at the time of exchange/sale
- It is the responsibility of the current equipment owner to ensure the service and inspection record form is kept up to date



•

8.6 PRODUCT CONFIGURATION

- Smirthwaite will document and maintain a record of the original product configuration at the time of first sale
- Smirthwaite will not be held responsible for any subsequent changes to this configuration unless authorized to do so in writing by Smirthwaite Ltd
- It is the equipment owner's responsibility to maintain their own records of changes to the equipment configuration and to be able to provide such records to subsequent owners to maintain traceability
- We recommend an inspection/service by a Smirthwaite Service Engineer (or Smirthwaite trained engineer) whenever a significant change is made to product configuration to ensure the product is safe to use. **If in any doubt, ALWAYS seek ADVICE**

	IMPORTANT REMINDER!
	DO NOT fit parts or accessories of other manufacturers to this product unless authorized to do so in writing by Smirthwaite Ltd.
	Any servicing or repairs required must be carried out by Smirthwaite Ltd (or a Smirthwaite trained engineer).
	If you believe this product to be faulty – DO NOT USE – Contact Smirthwaite Ltd on T: +44 (0) 1626 835552

If in any doubt, ALWAYS seek ADVICE.

8.7 EUROPEAN DIRECTIVE 2002/96/EC - WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT (WEEE)

The Waste Electrical and Electronic Equipment Directive (WEEE Directive) came into force in January 2007. The Directive aims to both reduce the amount of electrical and electronic equipment being produced and to encourage everyone to reuse, recycle and recover it.

As the responsible manufacturer, Smirthwaite operate a 'take back' scheme wherein we accept electrical and electronic components returned to our factory at the end of their product



service life.

Please note, we are obligated only to accept the return of electronic/electrical components for responsible recycling, not complete products or assemblies.

Please contact Smirthwaite Ltd on T: +44 (0) 1626 835552 should you require any further information.

8.0 WARRANTY & SERVICE

Smirthwaite Ltd warrants the products detailed on your order to be free from defects in materials and workmanship for a period of 2 years from date of delivery. If a fault develops during the period, please call Customer Services by email or telephone 01626 835552 who will advise you on the best course of action. Possible action may be for us to arrange to send out one of our Service Engineers, or have the goods returned to us. Should a repair not be possible within the guarantee period we will replace the product for new or nearest equivalent product. In the unlikely event that we cannot repair or exchange we will refund in full.

This warranty is for the UK only. The warranty excludes faults due to accident, neglect, misuse, not following the Instructions and normal wear and tear. This warranty is in addition to your legal rights. Goods will only be collected from the original delivery address.

A charge may be made where the goods cannot be repaired under the terms of the warranty. You will be advised before this is made.

T: +44 (0)1626 835552 E: info@smirthwaite.co.uk

10.0 CONTINUOUS IMPROVEMENT

Smirthwaite Ltd are committed to continuous improvement to their product range. Should you have any suggestions or comments please send them to our product design department at: info@smirthwaite.co.uk

Smirthwaite Ltd reserve the right to change the specification or material without prior notice.

For catalogues, help and further information on our products please contact us at:

Smirthwaite Ltd
16 Wentworth Road
Heathfield
Devon
TQ12 6TL

T: +44 (0)1626 835552

F: +44 (0)1626 835428



3.WHIRL ANGLE FOOTREST

INSTRUCTIONS FOR USE

Whirl Multi-Angle Footrest Code 8351-8354



CONTENTS	PAGE
1.0 INTRODUCTION	33
2.0 ILLUSTRATION OF YOUR FOOTREST	34
3.0 FOR YOUR SAFETY	35
4.0 UNPACKING AND ASSEMBLING YOUR FOOTREST	36
5.0 SETTING UP AND ADJUSTING YOUR FOOTREST	37
6.0 CARE & MAINTENANCE	39
7.0 WARRANTY & SERVICE	42
8.0 CONTINUOUS IMPROVEMENT	42
9.0 SERVICE INSPECTION	43

1.0 INTRODUCTION

Thank you for choosing the Smirthwaite Whirl Multi-Angle Footrest.

The Whirl Angle Footrest has been specifically designed to allow for independent multi-angle adjustment of the feet and ankles.



IMPORTANT!

These instructions should be read by all therapists and carers using the equipment and should be retained for future reference. The product should always be used under adult supervision.


Any incorrect use of the product and failure to follow the instructions may put the user at risk or impede the function. If you have any queries using this product or wish for further copies, please do not hesitate to contact Customer Service department on T: +44 (0) 1626 835552.

This product has been designed and manufactured specifically and solely for use by clients with special needs.



2.0 ILLUSTRATION OF YOUR FOOTREST

- A – RH Footplate
- B – LH Footplate
- C – Main Extrusion
- D – Footrest Bracket Assembly



3.0 FOR YOUR SAFETY

	STOP!
Please read these instructions CAREFULLY and THOROUGHLY	

- The user should **NOT** be left unattended whilst using the equipment. Always ensure a responsible therapist or carer is in attendance.
- The carer should be familiar with the methods of adjustment and have completed all adjustments appropriately to meet the needs of the child before transferring the child onto the equipment.
- Regular maintenance checks and cleaning are essential for the safe use of this equipment (see care and maintenance section).
- If any part is loose, damaged or functioning incorrectly, **DO NOT** use until rectified.
- **DO NOT** fit parts or accessories of other manufacturers to this product unless authorised to do so in writing by Smirthwaite Ltd. Failure to follow these instructions will not only invalidate the guarantee but could make the equipment dangerous to use. Smirthwaite Ltd will not accept liability for any injury or damage incurred through such malpractices. Any repairs required must be carried out by Smirthwaite Ltd authorised personnel.
- If you believe the equipment or any fitted accessory to be faulty at any time, **DO NOT USE** – contact Smirthwaite by telephone on +44 (0)1626 835552



4.0 UNPACKING AND ASSEMBLING YOUR FOOTREST






STOP!

If in any doubt, **ALWAYS** seek **ADVICE**.

Always turn hand wheels, levers and screws clockwise to tighten or anti-clockwise to loosen.

- When delivered, the footrest will be supplied assembled with the following fasteners ready to attach to your relevant product:
 - *2x T-Slot Nuts 8 St M8 (S02958)*
 - *2x Screws M8x14 Socket Button Head BZP (S02531)*



- To assemble, firstly fit t-slot nuts into the slots either side of the main extrusion.
- Secondly, fit the footrest assembly over the extrusion by sliding it on from the bottom.
- Finally, secure in position with the screws.
- If any part is loose, damaged or functioning incorrectly – do not use until rectified.



STOP!

If you believe this product to be faulty – **DO NOT USE** – Contact Smirthwaite Ltd on T: +44 (0) 1626 835552




5.0 SETTING UP AND ADJUSTING YOUR FOOTREST


	<p>STOP! If in any doubt, ALWAYS seek ADVICE</p> <p>Always turn handwheels and screws clockwise to tighten or anti-clockwise to loosen. The footrest must be fully adjusted by a therapist or trained representative before use.</p>
	<p>IMPORTANT! It is recommended that the footrest should be used for a few days prior to any operation so that the patient can become familiar with it.</p> <p>It should then be checked in hospital after the operation, and set up again if necessary.</p> <p>The footrest can normally be used from an early postoperative stage.</p>

5.1 ADJUSTMENTS

- A – RH Angle Adjuster
- B – LH Angle Adjuster
- C – RH Pitch Adjuster (20mm Increments)
- D – LH Pitch Adjuster (20mm Increments)
- E – RH Tilt Adjuster (20mm Increments)
- F – LH Tilt Adjuster (20mm Increments)



G – Height Adjustment Screws (Both Sides)



6.0 CARE AND MAINTENANCE



IMPORTANT!

Cleaning is recommended on a regular basis

- Clean footrest with a damp cloth and mild detergent.
- Stubborn marks on the woodwork should be cleaned by using a soft brush.
- Do not soak or immerse the footrest in water.
- DO NOT use bleach, solvents, abrasives, synthetic detergents, wax polishes, antibacterial sprays or wipes.
- For further information please refer to MHRA or your NHS cleaning guidelines.

6.1 DAILY CHECKS

- Check the footrest to ensure no parts are loose or damaged.
- Check all parts for signs of wear and tear.
- Check that the hand wheels tighten correctly, and are not cross-threaded or damaged.



STOP!

The user should not be using the equipment while the checks are carried out.

6.2 SERVICE INTERVAL

The footrest should be serviced every 3 years. Servicing must only be undertaken by a Smirthwaite service engineer, or by a Smirthwaite trained representative.

6.3 NOMINAL SERVICE LIFE

Your product has a nominal service life of 12 years, during which full post-sales support will be available with regard to spares and servicing.

Product service life has been determined based upon the design complexity of product, and the anticipated exposure to normal use.

Good practice dictates all Smirthwaite products have been designed and manufactured to high levels of safety and quality, and will meet requirements of normal use when maintained in line with our servicing recommendations.





STOP!

If the product has been out of use for an extended period of time (6 months or more) it should always be serviced prior to being reissued.

If the product has been subjected to 'heavy' or 'constant' use, the service should be reduced to half the recommended period.

Constant and/or heavy use is considered to be:

- Daily use above 7 hours duration
- Weekly use above 5 days duration
- Monthly use above 10 months per year
- Use by a client who is at 90% to 100% of the maximum weight limit of the product. The maximum weight limit must **NEVER** be exceeded
- Use by a client who is extremely active, either voluntarily or involuntarily

6.4 EXTENDING NOMINAL SERVICE LIFE

At Smirthwaite we are proud to produce products that have a reputation for quality and durability.

We believe our products have the potential to provide benefits to our clients beyond the nominal service life documented above.

We will continue to provide full support beyond the nominal service life provided the following conditions are met:

1. A full service schedule has been maintained.
2. A full service and inspection is undertaken at the end of the nominal service life period
3. The product is subsequently serviced annually (or biannually if under 'heavy/constant' use conditions)
4. Smirthwaite reserve the right to limit support where parts/components are no longer available

6.5 DOCUMENTATION/RECORDS

- It is the responsibility of the current equipment owner to ensure the 'Instructions for Use' manual and any further manuals for accessories fitted to the equipment are handed over to the new owner at the time of exchange/sale



- It is the responsibility of the current equipment owner to ensure the service and inspection record form is kept up to date

6.6 PRODUCT CONFIGURATION

- Smirthwaite will document and maintain a record of the original product configuration at the time of first sale
- Smirthwaite will not be held responsible for any subsequent changes to this configuration unless authorized to do so in writing by Smirthwaite Ltd
- It is the equipment owner's responsibility to maintain their own records of changes to the equipment configuration and to be able to provide such records to subsequent owners to maintain traceability
- We recommend an inspection/service by a Smirthwaite Service Engineer (or Smirthwaite trained engineer) whenever a significant change is made to product configuration to ensure the product is safe to use. **If in any doubt, ALWAYS seek ADVICE**

	<p>IMPORTANT REMINDER!</p> <p>DO NOT fit parts or accessories of other manufacturers to this product unless authorized to do so in writing by Smirthwaite Ltd.</p> <p>Any servicing or repairs required must be carried out by Smirthwaite Ltd (or a Smirthwaite trained engineer).</p> <p>If you believe this product to be faulty – DO NOT USE – Contact Smirthwaite Ltd on T: +44 (0) 1626 835552</p> <p>If in any doubt, ALWAYS seek ADVICE.</p>
---	--



7.0 WARRANTY & SERVICE

Smirthwaite Ltd warrants the products detailed on your order to be free from defects in materials and workmanship for a period of 2 years from date of delivery. If a fault develops during the period, please call Customer Services by email or telephone 01626 835552 who will advise you on the best course of action. Possible action may be for us to arrange to send out one of our Service Engineers, or have the goods returned to us. Should a repair not be possible within the guarantee period we will replace the product for new or nearest equivalent product. In the unlikely event that we cannot repair or exchange we will refund in full.

This warranty is for the UK only. The warranty excludes faults due to accident, neglect, misuse, not following the Instructions and normal wear and tear. This warranty is in addition to your legal rights. Goods will only be collected from the original delivery address.

A charge may be made where the goods cannot be repaired under the terms of the warranty. You will be advised before this is made.

T: +44 (0)1626 835552 E: info@smirthwaite.co.uk

8.0 CONTINUOUS IMPROVEMENT

Smirthwaite Ltd are committed to continuous improvement to their product range. Should you have any suggestions or comments please send them to our product design department at: info@smirthwaite.co.uk

Smirthwaite Ltd reserve the right to change the specification or material without prior notice.

For catalogues, help and further information on our products please contact us at:

Smirthwaite Ltd
16 Wentworth Road
Heathfield,
Devon.

TQ12 6TL

T: +44 (0)1626 835552

F: +44 (0)1626 835428

E: info@smirthwaite.co.uk

W: www.smirthwaite.co.uk



9.0 SERVICE INSPECTION

9.1 Product Information

Model	
Size	
Date of Manufacture	
Serial Number	
Final Inspection	

9.2 Service & inspection record form:

Date	Procedure	Service Personnel









Smirthwaite, 16 Wentworth Road, Heathfield, Newton Abbot, Devon.TQ12 6TL
T: +44 (0) 1626 835552 F: +44 (0) 1626 835428 E: info@smirthwaite.co.uk www.smirthwaite.co.uk

