



NEW SEATING 2026

# Which Asiento seat is right?

Three seat cushions, each matched to how the person sits. Start with what you see in the assessment.

Slides forward, or sits in posterior pelvic tilt?



**Ramped**

Needs lower-limb support to hold a stable midline?



**Contoured**

An active sitter who moves to reach and engage?



**Flat**

See each cushion in detail overleaf →

**smirthwaite**  
Together we  
transform lives

# Choosing the right seat cushion

Every Asiento is built around the pelvis. Choose the seat that gives the right foundation, then add only the support the person needs.

## Ramped

Holds the pelvis

Secures the pelvis in a neutral, aligned position and prevents forward sliding.

**Consider when** the user slides forward or loses pelvic position during the day, posterior pelvic tilt is present, or posture drops after 20–30 minutes.

**Pair with** sacral pad and 4-point pelvic belt.

## Contoured

Our standard seat

Increased immersion and lateral support to the lower limbs for a stable, midline position.

**Consider when** the legs drift out of alignment or rotate, midline sitting is hard to hold, or the therapist wants a shaped base for stability.

**Pair with** hip pads and abduction/adduction pads.

## Flat

Freedom to move

A neutral base for active users, supporting functional movement and greater independence.

Consider when the user moves to reach, play and engage, and the priority is freedom rather than containment.

Pair with only what is needed: hip, knee or lateral pads can be added to build support.

Cushion choice should always be confirmed through a clinical seating assessment. All three cushions work with the full Asiento accessory range.

## Not sure which seat?

Book an Asiento assessment with our team.

+44 (0)1626 835 552  
[info@smirthwaite.co.uk](mailto:info@smirthwaite.co.uk)  
[smirthwaite.com](https://smirthwaite.com)