

★asiento
★indiGO
★hardrock
★brookfield
★mobi

introducing our

NEW

plus range



This bespoke range offers a selection of products to meet the unique needs of plus sized individuals in both school and home environments

supporting children &
their families for life

smirthwaite



introducing our plus range

Recognising the increasing need for bespoke plus solutions, Smirthwaite is proud to introduce a range of products designed to enhance plus sized individuals' quality of life, making clear a difference by contributing to happier and more fulfilling lives.

From plus sized postural support chairs to heavy duty changing benches with larger weight and height requirements. Our plus range offers a selection of solutions to meet the unique needs of plus sized individuals in both school and home environments.

Individuals with larger body sizes (in height and/or weight) require specialised care and suitable equipment, posing challenges for occupational therapists, care facilities and providers. Limited resources affect timely care delivery, schooling and care home admissions, and hospital discharges.

Developed alongside input from clinical professionals, Smirthwaite's Plus range features a selection of plus products which are meticulously designed to enhance the quality of life for both users and carers, promoting equitable access to care and support.



★ asiento

3 designed specifically for plus sized children, teenagers and adults with moderate to complex needs

key features

- ★ **Durable** - accessories attach to a strong central column.
- ★ Optional **dynamic back** adjusts pressure for desired movement, aiding users needing sensory & physical feedback.
- ★ **Closer contact** - open chassis design & pushbike handles enable carers to get close to the backrest, beneficial for upper-limb therapy, care administration & calming users after an epileptic fit or during anxiety.



key benefits

- ★ Created to achieve improved postural support and symmetry, the Asiento delivers many benefits to the user and those around them, including:
 - enhanced eye contact and communication
 - improved interaction, breathing and eating
 - greater participation in play and activities.
- ★ **Levitex** technology is available as an optional extra.



6 designed for plus sized individuals with mild to moderate postural needs, who want to achieve enhanced independence

key features

- ★ **Adjustable lumbar support** for comfort.
- ★ **Flip away armrests** for easy self-transferring.
- ★ **Contoured seat** for proper pelvic positioning.
- ★ **Adjustable thoracic support** for customised comfort & posture management.



key benefits

- ★ This smart and stylish chair is designed for teenagers for both school and home use, providing **comfort, postural support, and integration** into any classroom setting.
- ★ It **promotes independence** with user-operated brakes and customisation options in black or blue, **size 6** for larger individuals.
- ★ Available with an **electric hi-lo base** for easy height adjustments, catering to different postural needs.
- ★ **Levitex** technology is available as an optional extra.



5

a robust & customisable chair for plus sized individuals available with a high or low back

key features

- ★ **Moulded back** to promote central alignment.
- ★ **Contoured seat** for proper pelvic positioning.
- ★ **Buckle pelvic strap** provides firm contoured support to stabilise hip position.
- ★ **Height adjustable arm rests** for customised comfort & posture.



key benefits

- ★ **Increased stability** – robust chair from 2 – young adults who require moderate levels of postural support.
- ★ **Greater adjustability** - seat height and depth adjustability. As an additional feature, the back can also be adjusted in height. This allows even greater flexibility & enables the chair to grow with your child.
- ★ **Perfect for customisation** - each child is unique, which is why the Brookfield has been carefully designed to provide the plus individual with a high level of customisation. By using the chair and the accessories together we can offer a highly tailored plus solution.
- ★ **Levitex technology** is available as an optional extra.



* **hardrock**

5 an extremely durable chair for children & adults who are prone to rock

key features

- ★ **Tray** provides trunk support & can be used for activities, playing & feeding.
- ★ **Anti thrust cushions** provides secure postural management.
- ★ **4-point pelvic strap** provides firm support to stabilise hip position.
- ★ **Mobile base & push handle** enables individual to be easily moved whilst seated.



key benefits

- ★ **Robust & durable** – specifically designed for children to young adults who are prone to rock repetitively. It provides moderate levels of postural support and is suitable for both school and home environments.
- ★ **Tailored to child's needs** - it can be a challenge to engage children with complex learning difficulties in educational & everyday activities. The Hardrock tackles this issue and offers a positive solution. It provides a safe environment and one that is conducive to learning.
- ★ **Posture matters** - uses an anti-thrust system and anterior pad to provide secure postural management within clearly defined boundaries. Thoracic supports and stabilisers are available for children who exhibit sideways motion.
- ★ **Levitex technology** is available as an optional extra.

5

150
kg

167 - 188
cm



L

★ mobi heavy
duty changer

1

designed to meet the needs of plus sized users, includes an powerful motor and extra wide mattress

key features

- ★ Electronically operated.
- ★ Travel range of 370mm.
- ★ Weight capacity up to 500kg.
- ★ Mobile & easy to manoeuvre.
- ★ Stowaway side rails.
- ★ Dry room environments only.



key benefits

- ★ **An adult changing table** with larger weight requirements. Its powerful 12K Newton motor, coupled with the extra wide mattress platform, makes it the ideal solution.
- ★ **Does not require any wall mounting or floor fixing.** They are easy to manoeuvre which means the user can be transported from room to room with no trouble (door opening has to be 1100mm or more).
- ★ **Electronically operated**, which makes adjusting the height of the table simple, this is beneficial for both the carer and user.
- ★ **Levitex technology** is available as an optional extra.

1

500
kg

≡

L

why plus?

Recognising the increasing need for bespoke plus solutions, Smirthwaite is proud to introduce a range of products designed to enhance plus sized individuals' quality of life, making clear a difference by contributing to happier and more fulfilling lives. From plus sized postural support chairs to heavy duty changing benches with larger weight and height requirements, our plus range offers a selection of solutions to meet the unique needs of plus sized individuals in both school and home environments.

product specifications

	asiento	indiGO	brookfield	hardrock
Product code	1023	11116	6759	8205
Size		6	5	5
Age range	Plus sized users	Plus sized users	Plus sized users	Plus sized users
Max user weight (kg)	120kg	100kg	150kg	150kg
Seat depth (min-max mm)	365 - 600	440-520	450-600	450-600
Seat width (mm)	480	490	550	550
Seat height range (mm)	490 - 570	490-650	450-700	425-550
Back height (mm)	490 - 590	560-640	600-750	650 (high back)



mobi heavy duty changer

Powder Coated	Size	Table Dim		Travel Range	Safety Sides		Castor Centres		Nominal clearances	Max User Weight	Product Weight
		L mm	W mm	H mm	W mm	H mm	P mm	Q mm	G mm	kg	kg
3161	1	2000	960	550-920	1050	165	1500	800	150	500	100





smirthwaite

Together we
transform lives



Book your free assessment
& discover more about our
plus range



44 (0)1626 835552



info@smirthwaite.co.uk

Key



Product size



Home



Max. user weight



School



User height



Dry room



Levitex



supporting children &
their families for life

16 Wentworth Road, Heathfield,
Newton Abbot, Devon, UK, TQ12 6TL

**MADE IN
BRITAIN**