



Tables & Desks

seating • standing • toileting • bathing
changing • showering
desks • tables • sleeping • therapy
sensory • moving • handling

smirthwaite

Together we
transform lives

welcome

Since 1986, Smirthwaite has developed and retained its reputation as one of the leading designers and manufacturers of bespoke special needs equipment. In that time, we have dedicated ourselves to finding and providing effective solutions that offer significant benefits to babies, children and young adults in their daily lives. Working closely with children, parents, carers and therapists, we have developed an innovative and comprehensive range of seating, standing, tables, desks, toileting, bathing, changing, showering, moving & handling, sleeping, therapy and sensory products.



Desks & tables

Our wide range of children's desks & tables are designed to fit almost any type of home or school seating – including our own seats of course!

With most of our children's desks & tables you can tilt the surface, adjust the height and alter other factors until your table or desk is a perfect fit for your child's needs.

Our range of desks & tables have been designed to fit in to the school and home environment. The contemporary beauty of wood and the quality of the design enables our desks to stand the test of time.

We have products that are suitable for both stationary and mobile use and our Product Advisors can run through the features with you to make sure you end up with the best possible solution.

Innovation in our design

Special needs equipment is often overly complex, difficult to use and expensive to maintain. At Smirthwaite, we believe in well-made, easy to use products that stand the test of time. Through close consultation with leading therapists, parents and end users, we ensure our designs offer numerous clinical benefits whilst maintaining a customer-centric and innovative approach.

Warranty & servicing

All of our products come with a two year warranty as standard. An inspection, maintenance, service and repair facility for all Smirthwaite products is available to you, for complete peace of mind. This can be scheduled at a time that is convenient for you.

Call our Customer Services team today on +44 (0)1626 835552 or e-mail info@smirthwaite.co.uk to book your service, inspection or repair.

Easy to clean


All of our products are easy to clean. For instance, our wooden chairs have anti-microbial properties that ensure surfaces are kept hygienic. This helps to reduce any risk of cross infection and is especially useful in multi-user environments such as nurseries and schools.

All our upholstery is carefully chosen due to its clever material make up that contains antibacterial agents that help against the spread of MRSA, E coli and other bacteria.

contents

tables & desks

study desk	4
personal tilt desk	5
hi-lo desk	6
tables	9
colours	12





Perfect for school and home

Our range of desks have been designed to fit in to the school and home environment. The contemporary beauty of wood and the quality of the design enables our desks to stand the test of time.

Study Desk

The Study Desk has been cleverly designed to allow your child to take part in various activities. The desk also contains ample storage space, so books and toys can be hidden away after they have been used. The Study Desk can be used with a variety of chairs, particularly the Heathfield and the Brookfield.



study desk

quality of the design
enables our desks to
stand the test of time

Personal Tilt Desk

The Personal Tilt Desk is packed with additional features; the tilt functionality can transform the desk and allow your child to work on a variety of activities. The Personal Tilt Desk features a cut out and is highly adjustable, allowing easy access for a variety of chairs, including wheel chairs.



personal tilt desk

Study Desk technical data

all dimensions in mm

Code	Size	Height range	Width	Depth	Max weight (kg)	Product weight (kg)
6967	1	345-520	450	600	80	19
6968	2	520-820	550	750	80	19

Personal Tilt Desk technical data

all dimensions in mm

Code	Size	Type	Height range	Depth	Width	Cut-out	Max weight (kg)	Product weight (kg)
6894	1	static	465-715	600	685	410	80	15
6895	1	mobile	540-790	600	685	410	80	15
6896	2	static	465-715	600	805	410	80	16
6897	2	mobile	540-790	600	805	410	80	16



hi-lo desk

Quick Height Adjusting Desks

Our range of desks have been designed to fit into the school environment with larger sized desk tops for dynamic learning.

Each of the desks have a modern metal mechanism with a turn handle to easily change the height of the table to the desired height, ideal when there are multiple people having use of the table. Each top has a sliding bottom extension to secure pens and papers from sliding off the desk when they are being used in a tilted position.

Our desks can be used with a variety of chairs and can be made to match the colour of our Juni chairs.

designed to fit into the school environment with larger sized desk tops for dynamic learning

Soft Close Tilting Desks


The hi-lo Deluxe Tilt Desk has a metal tilting mechanism and includes soft close hinges. Adjust the height quickly with a simple turn handle.

hi-lo deluxe tilt desk

Colours



Deluxe soft close hinges



Pen stop



buy online

Hi-lo fixed top

Shown with Juni size 3 on mobile skis



Desk pictured shown with optional coloured lime top. Standard colour is natural wood, unless otherwise requested on the order.

hi-lo fixed top


all dimensions in mm unless otherwise stated, price for standard natural wood top

Code	Height range	Width	Depth	Max weight (kg)	Product weight (kg)
2500	470 - 610	900	600	50	40
2501 (with castors)	545 - 685	900	600	50	40
2504	470 - 610	1200	600	50	45
2505 (with castors)	545 - 685	1200	600	50	45

hi-lo deluxe tilting top

all dimensions in mm unless otherwise stated, price for standard natural wood top

Code	Height range	Width	Depth	Max weight (kg)	Product weight (kg)
2502	470 - 610	900	600	50	40
2503 (with castors)	545 - 685	900	600	50	40
2506	470 - 610	1200	600	50	45
2507 (with castors)	545 - 685	1200	600	50	45



juni extras

Juni Table

The Juni table offers two key benefits; the legs are adjustable so your child can work or play at the correct height and the table folds, which allows for easy storage. Available in a variety of colours this is the perfect companion to your favourite Juni chair.



Juni table technical data
all dimensions in mm unless otherwise stated

Size	1	2
Height	310-450	400-600
Width	700	780
Depth	530	590
Max weight (kg)	60	60
Product weight (kg)	8	10



tables


range of tables and designs to suit most needs & requirements

Activity Table

Our Activity Table is an ideal height adjustable table for younger children at nursery, school and home. The steel legs can be adjusted in height to ensure good seating posture of your child and can easily be folded for storage.

Connectability

Our Connect range of rectangular and triangular tables have been designed to encourage good posture when your child is seated. All of the tables are adjustable in height and are extremely functional. They can be arranged into modular groups using rectangular and triangular configurations. Our Connect Tables meet a variety of different classroom needs to ensure an inclusive environment for children.



activitytable



 buy online

tables

connect
triangular table



connect
rectangular table



Connect Rectangular Table technical data

all dimensions in mm

Code	Size	Height range	Width	Depth	Max weight (kg)	Product weight (kg)
2311	1	350-450	800	585	80	12
2312	2	425-575	800	585	80	12
2313	3	525-700	800	585	80	12

Connect Triangular Table technical data

all dimensions in mm

Code	Size	Height range	Width	Depth	Max weight (kg)	Product weight (kg)
2361	1	350-450	1170	830	80	12
2362	2	425-575	1170	830	80	12
2363	3	525-700	1170	830	80	12

Activity Table technical data

all dimensions in mm

Code	Size	Height range	Width	Depth	Max weight (kg)	Product weight (kg)
2251	1	385-595	750	600	80	12
2252	2	555-765	750	600	80	12




seating standing toileting bathing changing showering
desks tables sleeping therapy sensory moving handling

colours

bright and colourful

Laminate colour range



The colours above are available for the hi lo desk and juni tables, apart from the beech laminate.

Beech laminate is available for the study desk, personal tilt desk, connect triangular table, connect rectangular table and the activity table.

Colours are a guide only and may vary due to limitations of the print process. Material samples are available upon request. Colour range may change without prior notice, we try to colour match where possible. Please contact our customer support team for more information.

smirthwaite

Together we
transform lives

always here to help

UK Contact Details:

telephone

+44 (0)1626 835552

email

info@smirthwaite.co.uk

website

www.smirthwaite.co.uk

Smirthwaite Ltd

16 Wentworth Road, Heathfield,
Newton Abbot, Devon TQ12 6TL

Smirthwaite reserves the right to modify product specifications at anytime without prior notification. All measurements are approximations and can vary depending on the wood and materials used. All details shown are correct at time of print. Smirthwaite® is a registered trade mark of Smirthwaite Limited.

