

Asiento S3


smirthwaite
Together we
transform lives

Positioning Pad Range

Product & Dimension Guide

The Asiento S3 seating system offers a modular range of Hip, Knee and Lateral support pads, each available in multiple sizes to provide precise, comfortable postural support tailored to the individual user.

- Hip Pad Range – Large & Small
- Knee Pad Range – Large, Medium & Small
- Lateral Pad Range – Large, Medium & Small



Overview

The Asiento S3 positioning pads mount directly onto the Asiento S3 seating and standing frame, providing configurable hip, knee and lateral trunk support.


All pads are manufactured to Smirthwaite's healthcare-grade standards, with sizes selected to match the user's individual measurements and clinical requirements.

Dimensions below are nominal (mm); refer to full technical drawings for manufacturing tolerances.



01 Hip Pad Range

The Hip Pad provides lateral pelvic support and is available in two sizes. Both sizes share the same 190mm width, with the Large pad offering greater depth (150–160mm) for more substantial support, and the Small pad a more compact 100–110mm depth. Overall thickness is 25mm.





Size	Width	Depth	Thicknes
Large	190 mm	150 – 160 mm	25 mm
Small	190 mm	100 – 110 mm	25 mm



02 Knee Pad Range



The Knee Pad Range provides medial/lateral knee support and is available in three sizes. Large, Medium and Small pads share a common 190mm height; width reduces from 180mm (Large) to 140mm (Medium and Small). Each size features a mounting cut-out sized to suit its bracket position. Overall thickness is 25mm.



Size	Width	Height	Notch (W x H)	Thickness
Large	180 mm	190 mm	40 x 30 mm	25 mm
Medium	140 mm	190 mm	35 x 110 mm	25 mm
Small	140 mm	190 mm	30 x 70 mm	25 mm

03 Lateral Pad Range

The Lateral Pad provides trunk/lateral body support and is available in three sizes, sharing a common 90mm depth and 25mm thickness. Width varies from 120mm (Small) to 200mm (Large) to suit the required level of trunk support.



Size	Length	Depth	Thickness
Large	200 mm	90 mm	25 mm
Medium	160 mm	90 mm	25 mm
Small	120 mm	90 mm	25 mm

Notes

- All dimensions are nominal and given in millimetres (mm). Do not scale – refer to full technical drawings for manufacturing tolerances.
- All pads are supplied in dark navy upholstery as standard; alternative colours available on request.
- For fitting guidance, ordering codes and pricing, please contact the Smirthwaite design or sales team.