


W LEG TROUGH OPTIONS



CENTRE CELL


The Centre Cell can be manufactured as an individual cell to allow for retention and increased stuffing.

FIRM CENTRE CELL

To give greater abduction the Centre Cell configuration can have solid foam, as can any cell, in the Three Cell option.

THREE CELL

To give greater flexibility there is the three cell configuration, which with an empty lower section offers a much flatter element and can take solid fill if required in any of the upper three cells.



TWIN WRAP

All centre cell and three cell configurations can have the wrap sewn in the centre to give independent wraps for each leg. The Twin Wraps are made from Dartex. This addresses clients who cross their legs within the W Leg Trough.

WRAP MATERIALS

The standard mesh Wrap can be substituted for a Dartex one, this used if there is any worry over pressure on the knees.

OUTER SUPPORT FRAME

To introduce greater support and framing we can give the W Leg Trough a full outer cradle.



THERMAL REGULATION COVERS

All sizes can be fitted with a distinctive blue second outer cover to address Thermal Regulation, not only will this help with cooling but will also re-heat if needed.

ASK FOR DETAILS FOR
THERMAL REGULATION
PRESSURE MANAGEMENT
T ROLL OPTIONS

