

samba package

a versatile seating system
with a unique midline
support, suitable for mild
to complex conditions

Tailored to additional needs

The Samba is a versatile chair; it can cater for a wide array of conditions (from mild to complex) and can have as much or as little posture support as required by your child.

Package

Samba chair with standard base or electric hi-lo chassis.

Comes complete with:

- 5-piece foam roll kit
- 4-point midline pelvic strap
- new, easy to use tray
- ferndale style sandals
- rigid swing away lateral supports
- samba style headrest
- footrest
- chassis



Deluxe Colour Range



Samba technical data

all dimensions in mm unless otherwise stated

Size	1	2	3	4
Code	8801	8802	8803	8804
Tilt in space (°)	12 - 25	12 - 25	12 - 25	12 - 25
Seat height	310 - 605	310 - 605	310 - 605	310 - 605
Seat depth	180 - 280	245 - 340	300 - 400	375 - 500
Seat width	300	340	340	390
Lower leg length to footrest	150 - 170	200 - 250	255 - 335	310 - 460
Back angle (°)	6 - 18	6 - 18	6 - 18	6 - 18
Back height	360 - 400	420 - 460	470 - 540	480 - 560
Armrest height ¹	220-320	220-320	220-320	220-320
Armrest angle (°)	39	39	39	39
Armrest height adjustment	90	90	9-	90
Overall chair height ²	670 - 1005	740 - 1060	790 - 1110	910 - 1240
Base length and width	650 x 510	650 x 510	650 x 510	650 x 510
Tray size	580 x 530	625 x 530	625 x 530	675 x 585
Max user weight (kg)	80	80	80	80
Product weight (kg)	31	34	35	38

1 Measured from the floor

2 Measured from the floor to the top edge of the chair back

+44 (0)1626 835552

info@smirthwaite.co.uk

smirthwaite.co.uk

Smirthwaite Ltd
16 Wentworth Road, Heathfield,
Newton Abbot, Devon TQ12 6TL

