

supro packages

a supine stander
with a difference &
building confidence

Three in one

The Supro provides the option of supine, prone and vertical standing in one frame. Its design gives an angle range of 90° allowing it to achieve a horizontal and vertical position.


Package 1 (Size 1) Package 2 (Size 2)

Supro standard width board comes complete with:

- upholstered board
- lateral and pelvic supports
- knee cups or colour co-ordinated knee capillars
- popliteal foams half or full thickness
- soft head support
- cocoon sandals
- plain tray
- grab rail



Deluxe Colour Range



Supro technical data

all dimensions in mm unless otherwise stated

Size	1	2	1	2
Upholstered board width	Standard		Wide	
Code	8256	8257	8258	8259
Thoracic support (top) to foot plate height	600-1070	700-1250	600-1070	700-1250
Pelvic support (top) to foot plate height	430-870	450-1010	430-870	450-1010
Knee block (top) to foot plate height	220-570	200-590	220-570	200-590
Foot plate to floor height (excluding sandals)	75	75	70	70
Max width (between thoracic supports)	330	330	430	430
Tray size (depth and width)	582 x 637	582 x 637	585 x 740	585 x 740
Angle [°]	0 to 90	0 to 90	0 to 90	0 to 90
Base frame	1140 x 570	1140 x 570	1130 x 610	1130 x 610
Max user weight (kg)	85	85	85	85
Product weight (kg)	55	56	57	59

Instructional video available online

+44 (0)1626 835552
info@smirthwaite.co.uk
www.smirthwaite.co.uk

Supporting children &
their families for life

