

Stayflex™ MD Series Chest Support

Offers firm trunk control with extra freedom of movement in the shoulders. Dual-zone construction features materials that optimally balance stretch and resistance. A great solution for users with lordosis and others needing strong sternal support.

Patented, dual-zone stretch construction prevents harness from shifting up toward the neckline, minimizing risk of injury.



MD (multi-directional) Pull Straps can be configured as front-pull or rear-pull at the time of the fitting to meet user needs.



Swivel Buckles on the lower straps and Side-Release Buckles on the top straps facilitate harness removal.



Optional zipper makes removal even easier. Available on both Standard and Narrow Width styles.

Part Numbers



Size	Standard Width with Zipper	Standard Width without Zipper	Narrow Width with Zipper	Narrow Width without Zipper
XS	-	SH351XS	-	-
S	SH350S	SH351S	-	-
M	SH350M	SH351M	SH360M	SH361M
L	SH350L	SH351L	SH360L	SH361L
XL	SH350XL	SH351XL	SH360XL	SH361XL

Add suffix -B5 or -B6 to end of part number for selected end-fitting; see options at right.

End-Fitting Options

Includes four end-fittings. Choose between high-strength plastic Flat-Mounts or Cinch-Mounts and add suffix -B5 or -B6 to the end of the part number.

Plastic Flat-Mount™



- B5

Plastic Cinch-Mount™



- B6

CAUTION: Whenever an anterior trunk support is used, a properly-adjusted pelvic support should also be worn to prevent sliding.

Visit us at bodypoint.com
206.405.4555 or 800.547.5716 | sales@bodypoint.com | Fax: 206.405.4556

© 2020 Bodypoint. All Rights Reserved. BMM327 2020.2

