

# Stayflex™ MD Series Chest Support

Offers firm trunk control with extra freedom of movement in the shoulders. Dual-zone construction features materials that optimally balance stretch and resistance. A great solution for users with lordosis and others needing strong sternal support.

Patented, dual-zone stretch construction prevents harness from shifting up toward the neckline, minimizing risk of injury.



MD (multi-directional) Pull Straps can be configured as front-pull or rear-pull at the time of the fitting to meet user needs.



Swivel Buckles on the lower straps and Side-Release Buckles on the top straps facilitate harness removal.



Optional zipper makes removal even easier. Available on both Standard and Narrow Width styles.



## Standard Width without Zipper

- ..... XS
- ..... S
- ..... M
- ..... L
- ..... XL

## End-Fitting Options

Includes high strength plastic flat-mounts.

### Plastic Flat-Mount™



**CAUTION:** Whenever an anterior trunk support is used, a properly-adjusted pelvic support should also be worn to prevent sliding.