

smirthwaite 130 mobile hoist

This great value mobile moving and handling hoist boasts an impressive 130kg lifting capacity, powered leg opening and market leading Linak® actuators and motors.

This user-friendly mobile hoist supports the transfer/lift of mobility impaired patients with minimal strain or risk to the caregiver, while providing complete safety, dignity and comfort to the individual.



With an impressive 130kg lifting capacity, powered leg opening and market-leading Linak® actuators and motors, this robust unit is suitable for everyday use.

It has all the essential functions required to transfer a child or young adult safely in both professional and homecare environments. As the hoist is mobile, a single unit can care for many patients.

Why not pair this hoist with one of our popular slings?

Product Code

3515

Features & Benefits

- Cost effective solution without compromising on quality.
- Powder-coated steel frame, offering durability and stability. Up to 130kg safe working load.
- User-friendly Linak control system making it easy to raise/lower the carry bar and to operate the powered leg opening/closing.
- Hoist safety features offer peace of mind including emergency lowering, manual lowering and an emergency stop button.
- Easy to position with large swivel castors, ergonomic handle and electrically operated leg opening and closing.
- 5 year warranty