

FSG130



Weights and dimensions

SWL	130 kg
Unit weight	22.7 kg
Track weight	6.8 kg
Side post & foot weight	9.1 kg

A (horizontal length)	2445 – 3387mm
B (horizontal travel)	1900 – 2850 mm
C (internal clearance)	2102 – 3052 mm
D (height range)	2141 – 2507 mm
E (foot length)	1300 mm

Details

Warranty	5 years
Compatible hoist	CP130
Material	Anodized aluminum

Standards

IEC 62366-1:2015 ISO
10535:2006
EN ISO 15223-1:2016

This product is CE marked and
complies with EU MDR 2017/745

Contact

16 Wentworth Road, Heathfield,
Newton Abbot, Devon, TQ12 6TL

info@smirthwaite.co.uk
+44 (0)1626 835552

