



# Toileting & Bathing

seating • standing • toileting • bathing  
changing • showering  
desks • tables • sleeping • therapy  
sensory • moving • handling



supporting children and  
their families **for life**

# welcome

Since 1986, Smirthwaite has developed and retained its reputation as one of the leading designers & manufacturers of special needs equipment. In that time, we have dedicated ourselves to finding and providing effective solutions that offer significant benefits to babies, children and young adults in their daily lives. Working closely with children, parents, carers and therapists, we have developed a comprehensive range of seating, standing, tables, desks, toileting, bathing, changing, showering, moving & handling, sleeping, therapy and sensory products.



## Toileting & bathing

When a child can bathe themselves in the bath chair and go to the toilet themselves they receive a boost to their confidence which helps them to fit in and develop socially. Our toileting and bath chair products, suitable for children from 12 months and upwards, can also make life easier for parents. By reducing frustration, bath time can become a fun and engaging experience for everyone.

The toileting chairs are freestanding and supportive meaning they fit well in nursery, school and home environments. If more privacy is required our potty chair comes with high sides. The chairs come in a range of colours allowing for a different colour to their normal chair. The two variations of bath chairs are both lightweight and compact and can also fit in a wide range of baths.

## Innovation in our design

Special needs equipment is often overly complex, difficult to use and expensive to maintain. At Smirthwaite, we believe in well-made, easy to use products that stand the test of time. Through close consultation with leading therapists, parents and end users, we ensure our designs offer numerous clinical benefits whilst maintaining a customer-centric and innovative approach.

## Warranty & servicing

All of our products come with a two year warranty as standard. An inspection, maintenance, service and repair facility for all Smirthwaite products is available to you, for complete peace of mind. This can be scheduled at a time that is convenient for you.

Call our Customer Services team today on **+44 (0)1626 835552** or e-mail **info@smirthwaite.co.uk** to book your service, inspection or repair.

## Easy to clean

All of our products are easy to clean. For instance, our wooden chairs have anti-microbial properties that ensure surfaces are kept hygienic. This helps to reduce any risk of cross infection and is especially useful in multi-user environments such as nurseries and schools. All our upholstery is carefully chosen due to its clever material make up that contains antibacterial agents that help against the spread of MRSA, E coli and other bacteria.

# contents

## toileting & bathing

bath corner chair	4
combi	5
potty chair	8
seahorse plus	10
colours	15



# bath corner chair

 buy online

### More support

Our Bath Corner Chair is designed for children who require a little extra help with sitting balance whilst bathing. This is achieved by its angled back and hand rail, both of which provide additional security.

### Adjustable

Like the Bath Chair, our Bath Corner Chair provides numerous levels of adjustment. The pommel and pelvic strap can be adjusted to suit each individual child's needs.

### Compact

Lightweight and easy to clean, the Bath Corner Chair has a perfectly compact design. This allows it to sit low in the bath enabling your child to fully extend their legs in complete comfort. Its size also enables it to fit a wide range of baths.



### accessories



#### Colour range

White ☐

The Bath Corner Chair comes complete with velcro pelvic strap, adjustable pommel, telescopic back rail and telescopic handrail as standard.

#### Telescopic handrail

CODE HR012  
Holding the handrail helps to maintain seating balance

#### Bath Corner Chair technical data

all dimensions in mm unless otherwise stated

Size	1	2
Code	7511	7512
Seat height	75	75
Seat width	340	340
Seat depth	280	330
Back height	350	460
Max user weight (kg)	30	35
Product weight (kg)	3	4

“The Bath Corner Chair gives my child extra support and helps me to bathe her.”

Christine, Mum, Reading

# combi

### Progressive toileting

Our Combi toileting chair is a great starter chair. Its design enables it to be used as a supportive chair whilst at the same time providing numerous features to assist your child with the everyday toileting needs.

### Providing solutions

The Combi toileting chair's design also incorporates a potty for younger children. This can easily be removed to provide a commode seat for children who have progressed their toilet training.

### Adjustability

Good levels of adjustability are offered with the Combi toileting chair so that a wide range of needs can be met. Seat height and depth can be adjusted to suit the needs of your child.

The Combi toileting chair comes complete with potty, plastic commode seat, upholstered seat and velcro pelvic strap. You can choose from either our fixed arm or height adjustable arm Combi toileting chair.



### accessories



#### Splash guard

CODE SG001 - SG002  
Assists in leg abduction and keeps toileting contained  
Requires Toggle handrail



#### Toggle handrail

CODE HR014 - HR017  
Allows child to hold and maintain balance and to sit with minimal support



#### Higher backrest in lieu of standard

CODE BK172 - BK175  
Provides extra back support for taller children







### Standard key features



Fixed or height adjustable  
arm rests



Plastic commode seat &amp; pot



Moulded Potty Seat



Adjustable seat height &amp; depth

## Combi technical data

all dimensions in mm unless otherwise stated

Size	1	2	3	4
Code	7551 <sup>1</sup> /7556 <sup>2</sup>	7552 <sup>1</sup> /7557 <sup>2</sup>	7553 <sup>1</sup> /7558 <sup>2</sup>	7554 <sup>1</sup> /7559 <sup>2</sup>
Seat height <sup>3</sup>	215-290	275-375	300-400	340-450
Seat width	290	325	365	435
Seat depth	205-290	245-355	310-385	350-440
Back height	325	370	400	500
Max user weight (kg)	35	45	55	65
Product weight (kg)	11	13	14	15

<sup>1</sup> Adjustable arm / <sup>2</sup> Fixed arm

<sup>3</sup> Seat height is to the top of upholstered seat (deduct 25mm for the height of the potty and deduct 40mm for the height of the commode seat)



## combi

	Size 1	Size 2	Size 3	Size 4
Lower leg accessories	Code	Code	Code	Code
skis	SK051	SK052	SK053	SK054
ski & half footrest	SK041	SK042	SK043	SK044
ski & full footrest	SK037	SK038	SK039	SK040
reversible footboard	FB023	FB024	FB025	FB026
adjustable footrest	FR010	FR011	FR012	FR013
mobile dolly base	DL012	DL013	DL014	DL015
footrest & straps¹	FP002	FP003	FP004	FP005
4 glisdomes complete with screws	GD001	GD001	GD001	GD001
Handrail accessories				
splash guard	SG001	SG002	SG002	SG002
toggle handrail	HR014	HR015	HR016	HR017
tray	TR006	TR073	TR074	TR005
Back accessories				
higher backrest in lieu of standard	BK175	BK172	BK173	BK174
adjustable height push handle	PH028	PH006	PH007	PH008
fixed height push handle	PH029	PH012	PH013	PH014
Cushions for adjustable arm Combi only				
25mm pelvic cushions (pair)	CS005	CS006	CS007	CS008
38mm pelvic cushions (pair)	CS014	CS015	CS016	CS017
50mm pelvic cushions (pair)	CS023	CS024	CS025	CS026
62mm pelvic cushions (pair)	CS032	CS033	CS034	CS035
75mm pelvic cushions (pair)	CS041	CS042	CS043	CS044
25+25mm pelvic cushions (pair)	CS050	CS051	CS052	CS053
38+25mm pelvic cushions (pair)	CS231	CS232	CS233	CS234
50+25mm pelvic cushions (pair)	CS077	CS078	CS079	CS080
62+25mm pelvic cushions (pair)	CS086	CS087	CS088	CS089
75x50mm wedge cushions (pair)	CS095	CS096	CS097	CS098
75x62mm wedge cushions (pair)	CS104	CS105	CS106	CS107
75x75mm wedge cushions (pair)	CS113	CS114	CS115	CS116
75x100mm wedge cushions (pair)	CS122	CS123	CS124	CS125

\*suitable for adjustable arm chairs only

1 Available on Mobile dolly base, Reversible footboard and Skis with full footrest

Combi comes complete with potty, plastic commode seat, upholstered seat and velcro pelvic strap.

Available with fixed arm or height adjustable arm chair.





# potty chair

## Fully supportive

Our Potty Chair has been designed specifically for young children who require assistance with sitting balance and need more reassurance when toileting. Its freestanding and highly supportive design means that it fits well in nursery, school and home environments. The sides of the Potty Chair provide your child with the privacy they sometimes need whilst toileting. Holding the handrail helps to maintain seating balance.

## Potty aperture

Many potties have a round aperture which can make toileting extremely difficult for children with slightly narrow hips. Our Potty Chair's aperture has been designed to ensure improved comfort and toileting performance.

## accessories



### Footboard

CODE FB012 – FB014  
Foot support that prevents chair movement. Factory fit at time of order



### Splash guard

CODE SG001 – SG002  
Assists in keeping legs apart and toileting contained  
Requires Toggle handrail

## Potty Chair technical data

all dimensions in mm unless otherwise stated

Size	1	2	3
Code	7525/7528 <sup>1</sup>	7526/7529 <sup>1</sup>	7527
Seat height	150-330	220-430	270-450
Seat width	310	380	475
Seat depth	330	380	430
User weight (kg)	35	45	55
Product weight (kg)	10	11	13

<sup>1</sup> Folding version (not available in size 3)

Colour Ranges  
See page 15 for colour options

freestanding and highly  
supportive design.  
The handrail helps to  
maintain sitting balance



The Potty Chair comes complete with potty, handrail and pelvic strap.





# seahorse plus

## Designed & Developed by NRS Healthcare

The Seahorse Plus inspires confidence and makes toileting and bathing a more manageable experience for both the carer and user.

Suitable for ambulant and non-ambulant users who require postural support, as well as users with mild to severe cognitive impairment or reduced motor function.

Fully adjustable for user security and comfort the Seahorse Plus can quickly and easily be positioned for use over a toilet or bath, as well as in the shower. It can also be used with a suitable hoist. Please contact the Customer Service team for advice on compatible hoists.

Available in three sizes: small, medium and large with a range of comfortable PU liners to accommodate children as they grow, from 5 years to adult.

A range of optional accessories are available, including; liners, harnesses and anterior supports.

Please note the Seahorse Plus is only available in blue.



a range of optional accessories have been designed around the user and carer to offer additional safety, security and maximum comfort



### Chest Harness:

Optional Poziform™ Chest Harness for users requiring more secure positioning when using the Seahorse Plus Toileting and Shower Chair.

**Small:** SH020

**Medium:** SH021

**Large:** SH022

**Extra Large:** SH023

# seahorse plus accessories

## Seahorse Plus Liners

Easy to fit, these soft PU foam liners have been designed for comfort and hygiene. The liners can be supplied in different thicknesses, allowing the chair to accommodate the user as they grow.

**Small Thick Liner:** SH013

**Small Thin Liner:** SH017

**Medium Thick Liner:** SH014

**Medium Thin Liner:** SH018

**Large Extra Thick Liner:** SH016

**Large Thick Liner:** SH015

**Large Thin Liner:** SH019



### Butterfly Harness:

Particularly suitable for women the optional Poziform™ Butterfly Harness offers more secure positioning when using the Seahorse Plus Toileting and Shower Chair.

**Small:** SH024

**Medium:** SH025

**Large:** SH026

**Extra Large:** SH027

## Product features



Tilt-in-Space: The innovative gas-assisted tilting mechanism enables 40° of rotation. To recline the seat simply press the red handle buttons then guide the user from a sitting position to fully reclined for safe, comfortable showering.





Seahorse Plus Flip-up Anterior Support  
(Replaces Armrests):

The Anterior Support helps with positioning and aids the user during the toileting process.

Small: SH044

Medium / Large: SH045

The Seahorse Plus has been created to support boys and girls as they grow. The chest harness is designed for boys and the butterfly harness for girls, for maximum comfort as they develop.

Small (SH001) Medium (SH002)  
Large (SH004)

Standard features	
Seahorse Plus Toileting & Shower Chair	
Height and depth adjustable headrest	
Height adjustable flip-up footrests	
Lap-strap	
Removableommel	
Removable armrests	
Removable commode pan	



seahorse plus  
dimensions



Small Seahorse Plus SH001							
Feature	Ref	Thick Liner		Thin Liner		No Liner	
		SH013		SH017			
Seat depth	SD	327mm	12 3/4"	355mm	14"	374mm	14 3/4"
Rear width of seat	SWR	258mm	10"	294mm	11 1/2"	332mm	13"
Front width of seat	SWF	294mm	11 1/2"	334mm	13 3/4"	365mm	14 3/4"
Height of back	BH	391mm	15 1/2"	391mm	15 1/2"	414mm	16 1/2"
Min headrest height	HH	456mm	18"	456mm	18"	475mm	18 1/2"
Max headrest height	HH	631mm	25"	631mm	25"	650mm	25 1/2"
Min height of seat	SH	540mm	21 1/2"	540mm	21 1/2"	521mm	20 1/2"
Max height of seat	SH	648mm	25 1/2"	648mm	25 1/2"	629mm	24 1/2"
Min height of footrests	FH	275mm	10 3/4"	275mm	10 3/4"	255mm	10"
Max height of footrests	FH	375mm	14 3/4"	375mm	14 3/4"	355mm	14"
Total width	OW	532mm	21"	532mm	21"	532mm	21"
Total depth (upright)	OD	807mm	31 3/4"	807mm	31 3/4"	807mm	31 3/4"
Max turning circle (titled)	TC	1209mm	47 1/2"	1209mm	47 1/2"	1209mm	47 1/2"
Min height over toilet	TH	461mm	18"	461mm	18"	461mm	18"
Max height over toilet	TH	569mm	22 1/2"	569mm	22 1/2"	569mm	22 1/2"
Min width over toilet	TW	426mm	16 3/4"	426mm	16 3/4"	426mm	16 3/4"
Armrest to seat	AH	179mm	7"	179mm	7"	199mm	8"
Seat angle	SA	0-40 °		0-40 °		0-40 °	
Max user weight		64kg	10st	64kg	10st	64kg	10st
Weight (including liners)		21kg	46.3lbs	20.2kg	44.5lbs	18.7kg	41.2lbs

Medium Seahorse Plus SH002							
No Liner		Thin Liner		Thin Liner		No Liner	
		SH014		SH018			
374mm	14 3/4"	425mm	16 3/4"	457mm	18 1/2"	457mm	18 1/2"
332mm	13"	303mm	12"	339mm	13 1/2"	339mm	13 1/2"
365mm	14 3/4"	353mm	14"	393mm	15 1/2"	393mm	15 1/2"
414mm	16 1/2"	490mm	19 1/2"	490mm	19 1/4"	490mm	19 1/4"
475mm	18 1/2"	562mm	22"	562mm	22"	562mm	22"
650mm	25 1/2"	737mm	29"	737mm	29"	737mm	29"
521mm	20 1/2"	521mm	20 1/2"	521mm	20 1/2"	521mm	20 1/2"
629mm	24 1/2"	625mm	24 1/2"	625mm	20 1/2"	625mm	20 1/2"
255mm	10"	325mm	12 3/4"	325mm	24 1/2"	325mm	24 1/2"
355mm	14"	425mm	16 3/4"	425mm	16 3/4"	425mm	16 3/4"
532mm	21"	586mm	23"	586mm	23"	586mm	23"
807mm	31 3/4"	860mm	33 3/4"	860mm	33 3/4"	860mm	33 3/4"
1209mm	47 1/2"	1413mm	55 1/4"	1413mm	55 3/4"	1413mm	55 3/4"
461mm	18"	442mm	17 1/2"	442mm	17 1/2"	442mm	17 1/2"
569mm	22 1/2"	551mm	21 3/4"	551mm	21 3/4"	551mm	21 3/4"
426mm	16 3/4"	214mm	18 3/4"	478mm	18 3/4"	478mm	18 3/4"
199mm	8"	179mm	8 1/2"	214mm	8 1/2"	214mm	8 1/2"
0-40 °		0-40 °		0-40 °		0-40 °	
64kg	10st	102kg	16st	102kg	16st	102kg	16st
18.7kg	41.2lbs	26.1kg	57.5lbs	24.7kg	54.5lbs	24.7kg	54.5lbs



# seahorse plus

## dimensions

Large Seahorse Plus SH003									
Feature	Ref	Extra Thick Liner SH016		Thick Liner		Thin Liner		No Liner	
Seat depth	SD	470mm	18 1/2"	455mm	18"	455mm	18"	455mm	19 1/4"
Rear width of seat	SWR	378mm	14 3/4"	411mm	16"	447mm	17 1/2"	447mm	19"
Front width of seat	SWF	420mm	16 1/2"	460mm	18"	500mm	19 3/4"	500mm	20 3/4"
Height of back	BH	549mm	21 1/2"	549mm	21 1/2"	549mm	21 1/2"	549mm	22 1/2"
Min headrest height	HH	584mm	24 1/2"	620mm	24 1/2"	620mm	24 1/2"	620mm	25"
Max headrest height	HH	620mm	31 1/2"	795mm	31 1/2"	795mm	31 1/2"	795mm	32"
Min height of seat	SH	795mm	20 1/2"	521mm	20 1/2"	521mm	20 1/2"	521mm	19 1/2"
Max height of seat	SH	521mm	24 3/4"	629mm	24 3/4"	629mm	24 3/4"	629mm	24"
Min height of footrests	FH	629mm	12 3/4"	325mm	12 3/4"	325mm	24 3/4"	325mm	12"
Max height of footrests	FH	325mm	17 3/4"	450mm	17 3/4"	450mm	12 3/4"	450mm	17"
Total width	OW	450mm	27"	692mm	27"	692mm	27"	692mm	27"
Total depth (upright)	OD	692mm	34 1/2"	877mm	34 1/2"	877mm	34 1/2"	877mm	34 1/2"
Max turning circle (titled)	TC	877mm	58 1/2"	1488mm	58 1/2"	1488mm	58 1/2"	1488mm	58 1/2"
Min height over toilet	TH	442mm	17 1/2"	442mm	17 1/2"	442mm	17 1/2"	442mm	17 1/2"
Max height over toilet	TH	550mm	21 3/4"	550mm	21 3/4"	550mm	21 3/4"	550mm	21 3/4"
Min width over toilet	TW	592mm	23 1/4"	592mm	23 1/4"	592mm	23 1/4"	592mm	23 1/4"
Armrest to seat	AH	214mm	8 1/2"	214mm	8 1/2"	234mm	8 1/2"	234mm	9"
Seat angle	SA								
Max user weight		152kg	24st	152kg	24st	152kg	24st	152kg	24st
Weight (including liners)		29kg	63.9lbs	28.3kg	62.4lbs	24.8kg	61.1lbs	24.8kg	54.7lbs

### Accessories

	Small	Medium	Large	Extra Large
Thick liners	SH013	SH014	SH014	SH016
Thin liners	SH017	SH018	SH018	
Chest harness	SH020	SH021	SH021	SH023
Chest harness (no buckle)	SH028	SH029	SH029	
Butterfly harness	SH024	SH025	SH025	SH027
Anterior supports	SH044	SH045	SH045	



# colours



## bright and colourful

### Moxxa range - brushed/soft finish smooth vinyl

raspberry	grey	blue
orange	red	black

Moxxa range can be used on most products

### Just Colour range - vinyl

ocean blue	rouge	gingersnap	marigold	eden
cobalt	candy pink	cool blue	cherry blossom	lilac
sea green	citrus green	greylfriar	jasmine white	

Just Colour range can be used on most products

Colours are a guide only and may vary due to limitations of the print process. Material samples are available upon request. Colour range may change without prior notice, we try to colour match where possible.



seating standing toileting bathing changing showering desks tables  
sleeping therapy sensory moving handling

a helping hand around the World



always here to help

UK Contact Details:

telephone

+44 (0)1626 835552

email

info@smirthwaite.co.uk

website

www.smirthwaite.co.uk

Smirthwaite Ltd

16 Wentworth Road, Heathfield,  
Newton Abbot, Devon TQ12 6TL



Smirthwaite reserves the right to modify product specifications at anytime without prior notification.  
All measurements are approximations and can vary depending on the wood and materials used.  
All details shown are correct at time of print. Smirthwaite® is a registered trade mark of Smirthwaite Limited.

