



Therapy & Stools

seating • standing • toileting • bathing
changing • showering
desks • tables • sleeping • therapy
sensory • moving • handling



supporting children and
their families **for life**

welcome

Since 1986, Smirthwaite has developed and retained its reputation as one of the leading designers & manufacturers of special needs equipment. In that time, we have dedicated ourselves to finding and providing effective solutions that offer significant benefits to babies, children and young adults in their daily lives. Working closely with children, parents, carers and therapists, we have developed a comprehensive range of seating, standing, tables, desks, toileting, bathing, changing, showering, moving & handling, sleeping, therapy and sensory products.



Innovation in our design

Special needs equipment is often overly complex, difficult to use and expensive to maintain. At Smirthwaite, we believe in well-made, easy to use products that stand the test of time. Through close consultation with leading therapists, parents and end users, we ensure our designs offer numerous clinical benefits whilst maintaining a customer-centric and innovative approach.

Easy to clean

All of our products are easy to clean. For instance, our wooden chairs have anti-microbial properties that ensure surfaces are kept hygienic. This helps to reduce any risk of cross infection and is especially useful in multi-user environments such as nurseries and schools.

All our upholstery is carefully chosen due to its clever material make up that contains antibacterial agents that help against the spread of MRSA, E coli and other bacteria.

Warranty & servicing

All of our products come with a two year warranty as standard. An inspection, maintenance, service and repair facility for all Smirthwaite products is available to you, for complete peace of mind. This can be scheduled at a time that is convenient for you.

Call our Customer Services team today on **+44 (0)1626 835552** or e-mail **info@smirthwaite.co.uk** to book your service, inspection or repair.

Therapy & stools

Therapy and conductive education plays a central role in helping children to keep active, manage symptoms and recover from illness or injury.

We can provide therapy equipment suitable for many different types of therapy, helping children to unlock their potential and maximise their development.

We specialise in Conductive Education but can tailor products to suit a wide variety of therapies.

So what is Conductive Education?

Conductive Education (CE) is a holistic system of learning (hence education) for children and adults with neurological motor disorders which affects their ability to control movement and which can also cause a wide range of developmental difficulties. As well as motor difficulties, a further implications of damage to the central nervous system most likely will affect their social, communication and problem solving skills as well as their confidence. Conductive Education was initially developed for the needs of children with cerebral palsy although nowadays it is used to help all age group children and adults with other neurological disorders such as Parkinsons, multiple sclerosis, strokes and head injuries.

Conductive Education is an active and open lifelong approach for many children and adults that helps them living with their disabilities without fighting against them. It is a method of learning their own ways of doing 'things' and how to facilitate and manage their difficulties in a socially inclusive environment and then at any further stage of their lives.


The impact of CE can be seen not only in the improvement of their motor skills and functions but also in the development of their individual personalities and when they feel more positive about their future to lead independent, fulfilled lives at home, school or work.

Conductive Education it is not a cure – it is a long term process that can transform lives of many who with support of their families and friends are willing to stay committed, determined and motivated. Its ultimate goal is to allow them achieving their full potential in everyday situations today as well as in the future.

contents

therapy & stools

footrests, steps & bars	4
therapy bench & bolster	5
wheely stools	6
accessories	7
colours	8





footrests, steps & bars buy online

Interlocking steps CODE 4188



Interlocking steps

Our 3 Interlocking Steps come complete with a single handrail that can be locked into position. All 3 steps have the same width (625mm) and depth (370mm). The lower step is at a height of 125mm, the middle step at 250mm and the platform at 375mm. Max user weight of 100kg and product weight of 20kg.

An extra handrail (CODE 4189) for the Interlocking steps is available to buy separately.

improving your child's balance

Parallel bars CODE 4292 - 4293

Our Parallel Bars are height and width adjustable and aim to improve your child's balance. They are extremely easy to use and adjust.



Interlocking footrest CODE 4185-4186

Our Interlocking Footrest is portable and can be used with a number of seating or therapy products. There are three interlocking footrests and they can be used together or separately. Available in 2 sizes.



Parallel Bars technical data

all dimensions in mm unless otherwise stated

Size	1	2
Code	4292	4293
Length	2000	2660
Height min	360	500
Height max	540	820
Width min	380	500
Width max	600	730
Max user weight (kg)	80	100
Product weight (kg)	15	25

Interlocking Footrest technical data

all dimensions in mm unless otherwise stated

Size	1	2
Code	4185	4186
Width	390	420
Depth	340	300
Height min	50	75
Height max	140	230
Heights individually	50 - 63 - 75	75 - 88 - 100
Max user weight (kg)	80	80
Product weight (kg)	5	5



therapy bench

buy online

“For our children, we look for products that will have great clinical benefits for the management of their posture.”

Anna, Teacher, Midlands

Our popular height and angle adjustable Therapy Bench together with accessories¹ is designed to assist in a wide range of children’s therapy programmes both at home and school.

Therapy Bench technical data

all dimensions in mm unless otherwise stated

Size	1	2	3	4	5	6	7	8
Code	5541	5542	5543	5544	5545	5546	5540	5547
Length	610	610	820	820	820	820	680	800
Width	250	300	300	380	325	380	280	360
Min height	195	195	280	280	405	405	170	270
Max height	295	295	455	455	705	705	270	400
Max user weight (kg)	100	100	100	100	100	100	100	100
Product weight (kg)	7	8	11	12	15	16	10	11



Colours



therapy bolster

buy online



Therapy bolster

CODE 5861 - 5866



Low end (single)

CODE LG014
(W 500mm x H 450mm)



High end (single)

CODE LG013
(W 500mm x H 700mm)

Colours



Our Therapy Bolster is extremely versatile and can be used just as a roll or built up to have a combination of High ends or Low ends.

Therapy Bolster technical data

all dimensions in mm unless otherwise stated

Size	1	2	3	4	5	6
Code	5861	5862	5863	5864	5865	5866
Length	600	1200	600	1200	600	1200
Height min	300	300	325	325	350	350
Height max	475	475	500	500	525	525
Diameter	200	200	250	250	300	300
Max user weight (kg)	100	100	100	100	100	100
Product weight (kg)	6	11	8	13	10	15

¹ Therapy Bench accessories, as shown in the image on the page opposite: Adjustable backrest with pelvic guides CODE **BK091** and Adjustable femoral kneeblocks (pair) CODE **KN002**

wheely stools

buy online



adjustable height wheely stool

Colours

- Red upholstery & yellow frame
- Green upholstery & blue frame
- Pink upholstery & blue frame



Therapist’s favourite

The Adjustable Height Stool is our most popular wheely stool. It has been designed specifically for working at lower levels with your child, student or client. The stool is designed for sitting astride the narrower side. Available in 3 colour options with a height range of 340-460mm. One size only.

Adjustable Height Wheely Stool technical data

all dimensions in mm

Code	Size	Colour	Width	Height min	Height max	Max user weight (kg)	Product weight (kg)
4181	1	Red & yellow	450	340	460	100	7
4182	1	Green & blue	450	340	460	100	7
4184	1	Pink & blue	450	340	460	100	7
4192	1	Red & yellow	500	340	460	150	8
4193	1	Green & blue	500	340	460	150	8
4194	1	Pink & blue	500	340	460	150	8

Round Wheely Stool technical data

all dimensions in mm

Code	Size	Colour	Height	Diameter	Base	Max user weight (kg)	Product weight (kg)
4178	1	Pacific blue	400	360	480 x 48	100	9
4179	1	Deep pink	400	360	480 x 48	100	9



round wheely stool

Colours




Space saving for therapists

Our space saving, Round Wheely Stool is a popular choice for therapists, teachers and parents when working at a low level with children. The stool comes in our natural wood finish and is available in one size only.



accessories




seating standing toileting bathing changing showering
desks tables sleeping therapy sensory moving handling



Suction hand grabs

CODE 4201-4203

Used to stabilise the hand to increase function in the opposing hand



Quick fix handrails

CODE 4205-4207

Will fit along, across or diagonally on any sized Slatted plinth



Raise blocks

CODE 4210-4214


A quick and simple method of raising the height of a plinth



Quick fix plastic mirror

CODE 4273 small (250 x 200mm)
CODE 4274 large (300 x 200mm)

Easy to attach and allows children to see themselves whilst using a plinth



Dycem non-slip mats

CODE 4703 red or blue (250 x 180mm)
CODE 4706 red or blue (350 x 250mm)
CODE 4709 red or blue (450 x 380mm)
CODE 4700 red or blue (round 190mm)




Handstick

CODE 4197 (single)
CODE 4198 (set of 5)

Can help to strengthen grip and increase dexterity


colours

Moxxa range - brushed/soft finish smooth vinyl



Moxxa range can be used on most products

Just Colour range - vinyl



Just Colour range can be used on most products

Colours are a guide only and may vary due to limitations of the print process. Material samples are available upon request. Colour range may change without prior notice, we try to colour match where possible.



always here to help

UK Contact Details:

telephone

+44 (0)1626 835552

email

info@smirthwaite.co.uk

website

www.smirthwaite.co.uk

Smirthwaite Ltd

16 Wentworth Road, Heathfield,
Newton Abbot, Devon TQ12 6TL

Smirthwaite reserves the right to modify product specifications at anytime without prior notification. All measurements are approximations and can vary depending on the wood and materials used. All details shown are correct at time of print. Smirthwaite® is a registered trade mark of Smirthwaite Limited.

