



# Toileting & Bathing

seating • standing • toileting • bathing  
changing • showering • desks • tables  
sleeping • therapy • sensory  
moving • handling • transportation



supporting children and  
their families **for life**

# welcome

Since 1986, Smirthwaite has developed and retained its reputation as one of the leading designers & manufacturers of special needs equipment. In that time, we have dedicated ourselves to finding and providing effective solutions that offer significant benefits to babies, children and young adults in their daily lives. Working closely with children, parents, carers and therapists, we have developed a comprehensive range of seating, standing, tables, desks, toileting, bathing, changing, showering, moving & handling, sleeping, therapy, transportation and sensory products.



## Toileting & bathing

When a child can bathe themselves in the bath chair and go to the toilet themselves they receive a boost to their confidence which helps them to fit in and develop socially.

Our toileting and bath chair products, suitable for children from 12 months and upwards, can also make life easier for parents. By reducing frustration, bath time can become a fun and engaging experience for everyone.

The toileting chairs are freestanding and supportive meaning they fit well in nursery, school and home environments. If more privacy is required our potty chair comes with high sides. The chairs come in a range of colours allowing for a different colour to their normal chair.

The two variations of bath chairs are both lightweight and compact and can also fit in a wide range of baths.

## Innovation in our design

Special needs equipment is often overly complex, difficult to use and expensive to maintain. At Smirthwaite, we believe in well-made, easy to use products that stand the test of time. Through close consultation with leading therapists, parents and end users, we ensure our designs offer numerous clinical benefits whilst maintaining a customer-centric and innovative approach.

## Warranty & servicing

All of our products come with a two year warranty as standard. An inspection, maintenance, service and repair facility for all Smirthwaite products is available to you, for complete peace of mind. This can be scheduled at a time that is convenient for you.

Call our Customer Services team today on **+44 (0)1626 835552** or e-mail **info@smirthwaite.co.uk** to book your service, inspection or repair.

## Easy to clean

All of our products are easy to clean. For instance, our wooden chairs have anti-microbial properties that ensure surfaces are kept hygienic. This helps to reduce any risk of cross infection and is especially useful in multi-user environments such as nurseries and schools.

All our upholstery is carefully chosen due to its clever material make up that contains antibacterial agents that help against the spread of MRSA, E coli and other bacteria.

# contents

## toileting & bathing

bath corner chair	4
combi	5
potty chair	8
seahorse plus	10
colours	15



# bath corner chair

 buy online

## More support

Our Bath Corner Chair is designed for children who require a little extra help with sitting balance whilst bathing. This is achieved by its angled back and hand rail, both of which provide additional security.

## Adjustable

Like the Bath Chair, our Bath Corner Chair provides numerous levels of adjustment. The pommel and pelvic strap can be adjusted to suit each individual child's needs.

## Compact

Lightweight and easy to clean, the Bath Corner Chair has a perfectly compact design. This allows it to sit low in the bath enabling your child to fully extend their legs in complete comfort. Its size also enables it to fit a wide range of baths.

## accessories



### Colour range

White

The Bath Corner Chair comes complete with velcro pelvic strap, adjustable pommel, telescopic back rail and telescopic handrail as standard.

### Telescopic handrail

CODE HR012

Holding the handrail helps to maintain seating balance

### Bath Corner Chair technical data

all dimensions in mm unless otherwise stated

Size	1	2
Code	7511	7512
Seat height	75	75
Seat width	340	340
Seat depth	280	330
Back height	350	460
Max user weight [kg]	30	35
Product weight [kg]	3	4

"The Bath Corner Chair gives my child extra support and helps me to bathe her."

Christine, Mum, Reading



# combi

## Progressive toileting

Our Combi toileting chair is a great starter chair. Its design enables it to be used as a supportive chair whilst at the same time providing numerous features to assist your child with the everyday toileting needs.

## Providing solutions

The Combi toileting chair's design also incorporates a potty for younger children. This can easily be removed to provide a commode seat for children who have progressed their toilet training.

## Adjustability

Good levels of adjustability are offered with the Combi toileting chair so that a wide range of needs can be met. Seat height and depth can be adjusted to suit the needs of your child.

The Combi toileting chair comes complete with potty, plastic commode seat, upholstered seat and velcro pelvic strap. You can choose from either our fixed arm or height adjustable arm Combi toileting chair.

## accessories



### Splash guard

CODE SG001 - SG002

Assists in leg abduction and keeps toileting contained  
Requires Toggle handrail



### Toggle handrail

CODE HR014 - HR017

Allows child to hold and maintain balance and to sit with minimal support



### Higher backrest in lieu of standard

CODE BK172 - BK175

Provides extra back support for taller children





 buy online


### Standard key features



Fixed or height adjustable arm rests



Plastic commode seat & pot



Moulded Potty Seat



Adjustable seat height & depth

# combi

	Size 1	Size 2	Size 3	Size 4
<b>Lower leg accessories</b>	<b>Code</b>	<b>Code</b>	<b>Code</b>	<b>Code</b>
skis	SK051	SK052	SK053	SK054
ski & half footrest	SK041	SK042	SK043	SK044
ski & full footrest	SK037	SK038	SK039	SK040
reversible footboard	FB023	FB024	FB025	FB026
adjustable footrest	FR010	FR011	FR012	FR013
mobile dolly base	DL012	DL013	DL014	DL015
footrest & straps <sup>1</sup>	FP002	FP003	FP004	FP005
4 glisdomes complete with screws	GD001	GD001	GD001	GD001
<b>Handrail accessories</b>				
splash guard	SG001	SG002	SG002	SG002
toggle handrail	HR014	HR015	HR016	HR017
tray	TR006	TR073	TR074	TR005
<b>Back accessories</b>				
higher backrest in lieu of standard	BK175	BK172	BK173	BK174
adjustable height push handle	PH028	PH006	PH007	PH008
fixed height push handle	PH029	PH012	PH013	PH014
<b>Cushions for adjustable arm Combi only</b>				
25mm pelvic cushions (pair)	CS005	CS006	CS007	CS008
38mm pelvic cushions (pair)	CS014	CS015	CS016	CS017
50mm pelvic cushions (pair)	CS023	CS024	CS025	CS026
62mm pelvic cushions (pair)	CS032	CS033	CS034	CS035
75mm pelvic cushions (pair)	CS041	CS042	CS043	CS044
25+25mm pelvic cushions (pair)	CS050	CS051	CS052	CS053
38+25mm pelvic cushions (pair)	CS231	CS232	CS233	CS234
50+25mm pelvic cushions (pair)	CS077	CS078	CS079	CS080
62+25mm pelvic cushions (pair)	CS086	CS087	CS088	CS089
75x50mm wedge cushions (pair)	CS095	CS096	CS097	CS098
75x62mm wedge cushions (pair)	CS104	CS105	CS106	CS107
75x75mm wedge cushions (pair)	CS113	CS114	CS115	CS116
75x100mm wedge cushions (pair)	CS122	CS123	CS124	CS125

\*suitable for adjustable arm chairs only

<sup>1</sup> Available on Mobile dolly base, Reversible footboard and Skis with full footrest

Combi comes complete with potty, plastic commode seat, upholstered seat and velcro pelvic strap.

Available with fixed arm or height adjustable arm chair.

### Combi technical data

all dimensions in mm unless otherwise stated

Size	1	2	3	4
Code	7551 <sup>1</sup> /7556 <sup>2</sup>	7552 <sup>1</sup> /7557 <sup>2</sup>	7553 <sup>1</sup> /7558 <sup>2</sup>	7554 <sup>1</sup> /7559 <sup>2</sup>
Seat height <sup>3</sup>	215-290	275-375	300-400	340-450
Seat width	290	325	365	435
Seat depth	205-290	245-355	310-385	350-440
Back height	325	370	400	500
Max user weight (kg)	35	45	55	65
Product weight (kg)	11	13	14	15

<sup>1</sup> Adjustable arm / <sup>2</sup> Fixed arm

<sup>3</sup> Seat height is to the top of upholstered seat (deduct 25mm for the height of the potty and deduct 40mm for the height of the commode seat)

 Colour Ranges  
See page 15 for colour options



Colour Ranges  
See page 15 for colour options



# potty chair

freestanding and highly supportive design. The handrail helps to maintain sitting balance

## Fully supportive

Our Potty Chair has been designed specifically for young children who require assistance with sitting balance and need more reassurance when toileting. Its freestanding and highly supportive design means that it fits well in nursery, school and home environments. The sides of the Potty Chair provide your child with the privacy they sometimes need whilst toileting. Holding the handrail helps to maintain seating balance.

## Potty aperture

Many potties have a round aperture which can make toileting extremely difficult for children with slightly narrow hips. Our Potty Chair's aperture has been designed to ensure improved comfort and toileting performance.


## accessories



### Footboard

CODE FB012 – FB014

Foot support that prevents chair movement. Factory fit at time of order



### Splash guard

CODE SG001 – SG002

Assists in keeping legs apart and toileting contained  
Requires Toggle handrail



## Potty Chair technical data

all dimensions in mm unless otherwise stated

Size	1	2	3
Code	7525/7528 <sup>1</sup>	7526/7529 <sup>1</sup>	7527
Seat height	150-330	220-430	270-450
Seat width	310	380	475
Seat depth	330	380	430
User weight (kg)	35	45	55
Product weight (kg)	10	11	13

<sup>1</sup> Folding version (not available in size 3)

Colour Ranges  
See page 15 for colour options

The Potty Chair comes complete with potty, handrail and pelvic strap.



# seahorse plus

## Designed & Developed by NRS Healthcare

The Seahorse Plus inspires confidence and makes toileting and bathing a more manageable experience for both the carer and user.

Suitable for ambulant and non-ambulant users who require postural support, as well as users with mild to severe cognitive impairment or reduced motor function.

Fully adjustable for user security and comfort the Seahorse Plus can quickly and easily be positioned for use over a toilet or bath, as well as in the shower. It can also be used with a suitable hoist. Please contact the Customer Service team for advice on compatible hoists.

Available in three sizes: small, medium and large with a range of comfortable PU liners to accommodate children as they grow, from 5 years to adult.

A range of optional accessories are available, including; liners, harnesses and anterior supports.

Please note the Seahorse Plus is only available in blue.



## Product features



**Tilt-in-Space:** The innovative gas-assisted tilting mechanism enables 40° of rotation. To recline the seat simply press the red handle buttons then guide the user from a sitting position to fully reclined for safe, comfortable showering.

a range of optional accessories have been designed around the user and carer to offer additional safety, security and maximum comfort



# seahorse plus accessories

## Seahorse Plus Liners

Easy to fit, these soft PU foam liners have been designed for comfort and hygiene. The liners can be supplied in different thicknesses, allowing the chair to accommodate the user as they grow.

**Small Thick Liner:** SH013

**Small Thin Liner:** SH017

**Medium Thick Liner:** SH014

**Medium Thin Liner:** SH018

**Large Extra Thick Liner:** SH016

**Large Thick Liner:** SH015

**Large Thin Liner:** SH019



### Chest Harness:

Optional Poziform™ Chest Harness for users requiring more secure positioning when using the Seahorse Plus Toileting and Shower Chair.

**Small:** SH020

**Medium:** SH021

**Large:** SH022

**Extra Large:** SH023



### Butterfly Harness:

Particularly suitable for women the optional Poziform™ Butterfly Harness offers more secure positioning when using the Seahorse Plus Toileting and Shower Chair.

**Small:** SH024

**Medium:** SH025

**Large:** SH026

**Extra Large:** SH027





**Seahorse Plus Flip-up Anterior Support  
(Replaces Armrests):**

The Anterior Support helps with positioning and aids the user during the toileting process.

**Small:** SH044

**Medium / Large:** SH045

The Seahorse Plus has been created to support boys and girls as they grow. The chest harness is designed for boys and the butterfly harness for girls, for maximum comfort as they develop.

**Small (SH001) Medium (SH002)  
Large (SH003)**

**Standard features**

- Seahorse Plus Toileting & Shower Chair
- Height and depth adjustable headrest
- Height adjustable flip-up footrests
- Lap-strap
- Removableommel
- Removable armrests



# seahorse plus dimensions

**Small Seahorse Plus SH001**

Features	Ref	Thick Liner	Thin Liner	No Liner
		SH013	SH017	
Seat depth	SD	327mm	355mm	374mm
Rear width of seat	SWR	258mm	294mm	332mm
Front width of seat	SWF	294mm	334mm	365mm
Height of back	BH	391mm	391mm	414mm
Min - max headrest height	HH	456mm - 631mm	456mm - 631mm	475mm - 650mm
Min - max height of seat	SH	540mm - 648mm	540mm - 648mm	521mm - 629mm
Min - max height of footrests	FH	275mm - 375mm	275mm - 375mm	255mm - 355mm
Total width	OW	532mm	532mm	532mm
Total depth (upright)	OD	807mm	807mm	807mm
Max turning circle (titled)	TC	1209mm	1209mm	1209mm
Min - max height over toilet	TH	461mm - 569mm	461mm - 569mm	461mm - 569mm
Min width over toilet	TW	426mm	426mm	426mm
Armrest to seat	AH	179mm	179mm	199mm
Seat angle	SA	0-40°	0-40°	0-40°
Max user weight		64kg	64kg	64kg
Weight (including liners)		21kg	20.2kg	18.7kg

**Medium Seahorse Plus SH002**

Features	Ref	Thick Liner	Thin Liner	No Liner
		SH014	SH018	
Seat depth	SD	425mm	457mm	470mm
Rear width of seat	SWR	303mm	339mm	378mm
Front width of seat	SWF	353mm	393mm	423mm
Height of back	BH	490mm	490mm	514mm
Min - max headrest height	HH	562mm - 737mm	562mm - 737mm	584mm - 759mm
Min - max height of seat	SH	521mm - 625mm	521mm - 625mm	501mm - 609mm
Min - max height of footrests	FH	325mm - 425mm	325mm - 425mm	305mm - 405mm
Total width	OW	586mm	586mm	586mm
Total depth (upright)	OD	860mm	860mm	860mm
Max turning circle (titled)	TC	1413mm	1413mm	1414mm
Min - max height over toilet	TH	442mm - 551mm	442mm - 551mm	442mm - 551mm
Min width over toilet	TW	478mm	478mm	478mm
Armrest to seat	AH	214mm	214mm	234mm
Seat angle	SA	0-40 °	0-40 °	0-40 °
Max user weight		102kg	102kg	102kg
Weight (including liners)		26.1kg	24.7kg	29kg



# seahorse plus dimensions

Medium Seahorse Plus SH003					
Features	Ref	Extra Thick Liner SH014	Thick Liner SH014	Thin Liner SH018	No Liner
Seat depth	SD	455mm	455mm	455mm	488mm
Rear width of seat	SWR	375mm	411mm	447mm	447mm
Front width of seat	SWF	420mm	460mm	500mm	530mm
Height of back	BH	549mm	549mm	549mm	572mm
Min - max headrest height	HH	620mm - 795mm	620mm - 795mm	620mm - 795mm	641mm - 816mm
Min - max height of seat	SH	521mm - 629mm	521mm - 629mm	521mm - 629mm	501mm - 609mm
Min - max height of footrests	FH	325mm - 450mm	325mm - 450mm	325mm - 450mm	305mm - 430mm
Total width	OW	692mm	692mm	692mm	692mm
Total depth (upright)	OD	877mm	877mm	877mm	877mm
Max turning circle (titled)	TC	1488mm	1488mm	1488mm	1488mm
Min - max height over toilet	TH	442mm - 550mm	442mm - 550mm	442mm - 550mm	442mm - 550mm
Min width over toilet	TW	592mm	592mm	592mm	592mm
Armrest to seat	AH	214mm	214mm	214mm	234mm
Seat angle	SA	0-40 °	0-40 °	0-40 °	0-40 °
Max user weight		152kg	152kg	152kg	152kg
Weight (including liners)		29kg	28.3kg	27.7kg	24.8kg

## Accessories

	Small	Medium	Large	Extra Large
Thick liners	SH013	SH014	SH014	SH016
Thin liners	SH017	SH018	SH018	
Chest harness	SH020	SH021	SH021	SH023
Chest harness (no buckle)	SH028	SH029	SH029	
Butterfly harness	SH024	SH025	SH025	SH027
Anterior supports	SH044	SH045	SH045	



# colours



bright and colourful

## Just Colour range - vinyl

ocean blue	rouge	gingersnap	marigold	eden
cobalt	candy pink	cool blue	cherry blossom	lilac
sea green	citrus green	greyfriar	jasmine white	

Just Colour range can be used on most products

Colours are a guide only and may vary due to limitations of the print process. Material samples are available upon request. Colour range may change without prior notice, we try to colour match where possible.





seating standing toileting bathing changing showering desks tables  
sleeping therapy sensory moving handling

a helping hand around the World



always here to help

UK Contact Details:

telephone

+44 (0)1626 835552

email

info@smirthwaite.co.uk

website

www.smirthwaite.co.uk

Smirthwaite Ltd

16 Wentworth Road, Heathfield,  
Newton Abbot, Devon TQ12 6TL



Smirthwaite reserves the right to modify product specifications at anytime without prior notification.

All measurements are approximations and can vary depending on the wood and materials used.

All details shown are correct at time of print. Smirthwaite® is a registered trade mark of Smirthwaite Limited.

