



Samba

a versatile paediatric seating
system with a unique midline
support, suitable for mild to
moderate conditions

Supporting children &
their families **for life**

why choose Samba?

Tailored to additional needs

The Samba is versatile chair, accommodating a range of needs from mild to complex, with adjustable posture support for children. It features a 4-point pelvic strap to ensure proper pelvic positioning, helping to prevent sacral sitting. The chair now also includes a new, user-friendly plastic tray.

Promoting postural control

To promote postural control for children needing support, the soft cross harness is recommended. It stabilizes the shoulders for spine control, unlike conventional laterals. For enhanced alignment, it can also be used with a 4-point pelvic strap.

A multi-height solution

It provides an impressive seat height range of 295mm. The seat height ranges from 310mm – 605mm, this makes it an ideal solution for children that study, play and live in a multi-height environment.



Understanding posture

The Samba is available in four sizes and has an unique approach to managing posture. It does this effectively by using the contours of the seat and backrest, as opposed to hard intervention.

The contouring can be tailored to the needs of the child by adding or taking away the foam inserts. The contours in the backrest have a subtle, but effective way of holding a child in a good mid-line position, without the intervention of using fixed laterals.

Do you have any bespoke paediatric standing requirements for children? Our friendly experts are here to help and offer **FREE assessments**.

get in touch

-  smirthwaite.co.uk
-  info@smirthwaite.co.uk
-  +44 (0)1626 835552

