



Transportation

Discover Smirthwaite's Paediatric Car Seat Range

Safety, comfort, and peace of mind on every journey

At Smirthwaite, we combine cutting-edge ergonomics, rigorous crash testing, and clinical insights to offer a comprehensive range of specialized car seats. Whether your child requires advanced physical support or specialised behavioral boundaries, our three flagship models ensure complete security on the road and total confidence for carers.

smirthwaite

SIMON

Advanced Postural Management & Stabilization



The SIMON is meticulously crafted to support children with complex physical needs, ensuring they maintain an optimal, comfortable position throughout the journey. Built on extensive research and crash testing, it seamlessly blends clinical support with innovative design.

- **First-Class Safety:** Fully guaranteed to the latest provisions of the stringent UNECE Regulation R129/03 standard.
- **Integrated Swivel & Tilt:** Features an integrated swivel base for easy transfers, alongside a built-in Tilt-in-Space system to optimize pressure management.
- **Grows with the Child:** Fully adjustable backrest height and seat depth to accommodate your child as they grow (suitable for user heights of 100cm to 150cm).
- **The "SOUND" System:** Innovative built-in headrest speakers allow children to connect devices and listen to calming music or audio during transit.

Best For: Children requiring high-level physical positioning, stabilisation, and custom adjustments.



Qai

Flexible Positioning & 360° Carer Convenience



The Qai is designed to maximise both child safety and carer ease-of-use. Authorised as a child restraint system for special needs (ECER129/04 S), it has undergone over 40 rigorous crash tests to meet the highest safety benchmarks.

- **360° Rotatable Base:** Unrestricted rotation provides effortless vehicle access, significantly reducing physical strain on parents and carers during transfers.
- **Flexible Facing Options:** Offers rearward-facing positioning for smaller children (61–105cm) and forward-facing options from 76cm up to 150cm.
- **Tailored, Adaptive Fit:** Includes rigid lateral supports, a comprehensive harness, and an abduction wedge, alongside removable inner cushions for smaller or slimmer children.
- **Adjustable Footrest:** Complete height and angle adjustability ensures proper, stable leg support.

Best For: Carers looking for seamless transfers and a highly adjustable, multi-facing growth seat.



Voyager

The Ultimate Solution for Behavioral & Postural Harmony

Manufactured by the industry-renowned In Car Safety Centre and distributed by Smirthwaite, the Voyager is widely recognised as the premier solution for children up to 36kg who require extra security, reassurance, and behavioral boundaries while travelling.

- **Behavioral Peace of Mind:** Engineered to provide a secure, comforting environment that prevents children from unsafely escaping their seat or harness during transit.
- **5-Point Safety Harness:** Equipped with a one-pull adjustment harness and a specialized chest clip to prevent shoulder straps from being slipped off.
- **Customisable Comfort:** The luxury cover accommodates extra foam rolls and side supports, with an optional recline bar for children requiring a deeper tilt.
- **Motability Scheme Approved:** Fully supported by the Motability Scheme, ensuring high accessibility for families across the UK.

Best For: Children with specialised behavioral needs, sensory processing differences, or autism who require secure boundaries.



Book a free assessment today

+44 (0) 1626 835552

info@smirthwaite.co.uk

smirthwaite.co.uk